SINCE 1966
MADE IN ITALY

NebulAIR

by **FLAEM**®

GB INSTRUCTIONS FOR USE MANUAL

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
E MANUAL DE INSTRUCCIONES

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Mod. NEBULAIR

Aerosol Therapy Device

Congratulations on your purchase and thank you for choosing our unit.

Our goal is to fully satisfy customers by offering them cutting-edge systems for the treatment of respiratory tract ailments. Carefully read these instructions and keep them for future reference. Use the unit only as described in this instruction manual. This is a medical device for home use and may only be **used with medication prescribed or recommended by a doctor**. Visit our Internet site www.flaem.it to view the whole range of Flaem products.

THE STANDARD COMPONENTS AND ACCESSORIES OF THE UNIT INCLUDE:

Ⓐ - Aerosoltherapy unit (main unit)

- A1 - On/off switch
- A2 - Air outlet port
- A3 - Air filter
- A4 - Nebuliser holder
- A5 - Carrying handle
- A6 - Cable compartment
- A7 - Power cord

Ⓑ - Connection tube (main unit / nebuliser)

Ⓒ - Accessories

- C1 - Adult mask
- C2 - Child mask
- C3 - Mouthpiece with valve
 - C3.1 Expiratory valve
- C4 - Non-invasive adult nasal prong
- C5 - Nebulisation manual control
- C6 - RF6 Plus nebuliser
 - C6.1 - Upper part
 - C6.2 - Complete nozzle
 - C6.3 - Lower part
 - C6.4 - Double action inspiratory valve
- C7 - Nebuliser and filter spares kit

⚠ IMPORTANT SAFEGUARDS


- Before using the device for the first time and periodically during the product's life, check the integrity of its structure and of the power cable, making sure that there is no damage; if it is damaged, do not insert the plug and take the product immediately to an authorised FLAEM technical assistance centre or to your local dealer.
- The average lifespan of accessories is 1 year; however, we recommend you replace the nebulizer every 6 months in intensive use (or earlier, if the nebulizer is blocked) to always ensure maximum therapeutic effect.
- Children and disabled persons should always use the unit under strict supervision of an adult who has read this manual.
- Several parts of the device are small enough to be swallowed by children; therefore keep the device out of reach of children.
- Do not use the tubes and the cords of the device for purposes other than those specified as they may cause strangulation hazard; utmost attention must be paid to children and persons with special needs; usually, such persons are not capable of assessing the hazard correctly.
- The unit should not be used in the presence of anaesthetic mixture inflammable with air, oxygen or nitrous oxide.
- Keep the cord away from heated surfaces.
- Keep the power supply cord away from animals (for example rodents) as they may damage the insulation of the power supply cord.

- Do not handle the unit with wet hands. Do not use the device in humid environments (e.g. while taking a bath or shower). Do not immerse the device in water; should this accidentally occur, unplug the device immediately. Do not remove or touch the device immersed in water before pulling out the plug. Take the device immediately to an authorized service center or to your local dealer.
- Only use the device in dust-free environments; otherwise, therapy may be compromised.
- The casing of this unit is not protected against the penetration of liquids. Do not wash the device under running water or submerge in water. Keep away from splashing water or other liquids.
- Do not expose the device to high temperatures.
- Do not put the device near heat sources, direct sunlight, or in hot environments.
- Do not obstruct or insert any objects in the filter or in its housing.
- Do not obstruct the air vents placed on both sides of the unit.
- During use, always place the unit on a firm surface clear of obstacles.
- Check that there is no material obstructing the air vents before each use.
- Do not place any objects inside the air vents.
- The average duration provided for compressor ranges is: F400: 400 hours, F700: 700 hours, F1000: 1000 hours, F2000: 2000 hours.
- Repairs must be performed by authorized personnel only. Any unauthorized repairs will void the warranty and may pose a safety hazard for the user.
WARNING: Do not modify this equipment without authorization of the manufacturer
- The Manufacturer, the Retailer and the Importer are responsible for the safety, reliability and performance of the device only if: a) the device is used in accordance with the operating instructions b) the electrical system of the premises where the device is used is in compliance with current laws.
- Interactions: The materials used that come into contact with medications have been tested with a wide range of medications. However, due to the variety and continuous evolution of medication, the possibility of a chemical interaction cannot be excluded. We recommend that once the medication has been opened you should use it as soon as possible and avoid prolonged contact of the medication with the nebulizer.
- The manufacturer should be contacted for reporting problems and/or unexpected events related to device operation.

INSTRUCTIONS FOR USE

Before each use, the nebuliser and accessories should be cleaned according to the instructions in the section “CLEANING, SANITIZATION, DISINFECTION AND STERILIZATION”. A personal use of the nebuliser and its accessories is recommended in order to avoid any risks of infection.

1. Insert the power cord (A7) into a power socket corresponding to the voltage of the device. This must be positioned so that it is not difficult to make the disconnection from the mains.
2. Open the nebuliser by turning the upper part (C6.1) anticlockwise.
3. Pour the drug prescribed by the physician into the lower part (C6.3). Close the nebuliser by turning the upper part (C6.1) clockwise.
4. Connect accessories as indicated in the “Connection diagram” on the cover.
5. Sit comfortably holding the nebuliser in your hand, place the mouthpiece onto your mouth or alternatively use the nose piece or mask. Should you opt for the

mask accessory, place it on your face as shown in the picture (with or without using the elasticated strap).

6. Start the unit by means of switch (A1) and breathe in and out deeply.

7. When the treatment is finished, turn off and unplug the unit.

WARNING: If after the therapy session moisture accumulates in the tube (B), remove the tube from the nebulizer and dry it using the device's compressor fan; this operation will prevent mold from forming inside the tube.

HOW TO USE RF6 PLUS NEBULISER WITH DOUBLE ACTION VALVE

It is professional, quick and suitable for administering all types of drugs, including more costly ones, also in patients with chronic diseases. The configuration of the internal lines of the RF6 Plus nebuliser ensures the ideal granulometry providing effective treatment right down to the lower airways.

Use the double action valve (C6.4) in position 1. In this case the inhalation therapy allows the maximum drug intake and therefore the maximum therapeutic effect, minimising the loss of medication.


To make inhalation therapy faster we recommend using the double action valve (C6.4) in position 2.


To administer the medication only through the nose, use the non-invasive nasal prong (C4).

USE OF NEBULISATION MANUAL CONTROL

To achieve continuous nebulisation action you should not use the manual nebulisation control (C5), especially in the case of children or persons with reduced physical, sensory, or mental capabilities. The manual nebulisation control is useful for limiting dispersion of the medication in the surrounding environment.

To start nebulizing close with a finger the hole of the nebulizer manual control (C5) and breathe in gently; we recommend to hold your breath for a moment so that the inhaled aerosol droplets can be deposited,

meanwhile, to disable nebulizing, remove your finger from the hole of the nebulizer manual control to avoid the waste of drug, optimizing its acquisition. Then exhale slowly.

NON-SPILL SYSTEM

The RF6 Plus nebuliser is also especially suitable for pediatric use, since the non-spill system prevents medication from spilling out even if the nebuliser is accidentally tipped over.

Efficacy tested up to 5 ml, the average dose recommended for normal use

CLEANING, SANITIZATION, DISINFECTION AND STERILIZATION

Turn off the unit before any cleaning and disconnect the power cord from the mains.

CLEANING THE DEVICE AND THE TUBE'S OUTER SURFACE

Only use a cloth dampened with antibacterial detergent (non-abrasive and free of solvents of any kind).

SANITIZATION AND DISINFECTION OF THE ACCESSORIES

We recommend a personal use of both the nebulizer and accessories to avoid any risk of infection or contagion. Open the nebulizer by turning the top (C6.1) counterclockwise, then remove the nozzle from the top (C6.1), and disassemble it as shown in the "Connection Diagram" section C6.

SANITIZATION

Before and after each use, sanitize the nebulizer and accessories according to one of the methods described below.

(Method A): The accessories that can be sanitized with the method A are (C1-C2-C3-C4-C5-C6.1-C6.2-C6.3-C6.4). Sanitize the accessories in warm (approximately 40° C) drinking water with mild dishwashing liquid (non abrasive) or in the dishwasher with hot cycle.

(Method B): The accessories that can be sanitized with the method B are (C1-C2-C3-C4-C5-C6.1-C6.2-C6.3-C6.4). Sanitize the accessories by immersion in a solution consisting of 50% water and 50% white vinegar.

Finally rinse thoroughly with warm (approximately 40° C) drinking water.

(Method C): The accessories that can be sanitized with the method C are (C3-C4-C5-C6.1-C6.2-C6.3-C6.4). Sanitize the accessories by boiling in water for 20 minutes; it is preferable the use of demineralized or distilled water to avoid scale deposits.

After sanitization shake the accessories vigorously and place them on a paper towel, or alternatively, dry them with a jet of hot air (e.g. using a hair dryer).

DISINFECTION

The accessories that can be disinfected are (C1-C2-C3-C4-C5-C6.1-C6.2-C6.3-C6.4).

The disinfection procedure described in this paragraph is to be performed before using the accessories and is effective only if all the prescriptions herein are met and all the components involved are previously sanitized.

The disinfectant to be used must be of electrolytic oxidizing chlorine type, (active ingredient: sodium hypochlorite), specific for disinfection and is available in all pharmacies.

Procedure:

- Fill a container, suitable to contain all the individual components to be disinfected, with a solution of drinking water and disinfectant, while respecting the proportions indicated on the packaging of the disinfectant itself.
- Completely immerse each individual component in the solution, taking care to avoid the formation of air bubbles in contact with the components. Leave the components immersed for the period of time indicated on the packaging of the disinfectant, and associated with the concentration chosen to prepare the solution.
- Remove the components now disinfected and rinse thoroughly with lukewarm drinking water.
- After disinfection shake the accessories vigorously and place them on a paper towel, or alternatively, dry them with a jet of hot air (e.g. using a hair dryer).
- Dispose of the solution as recommended by the disinfectant manufacturer.

STERILIZATION

The accessories that can be sterilized are (C3-C4-C5-C6.1-C6.2-C6.3-C6.4).

The sterilization procedure described in this paragraph is effective only if all the

prescriptions herein are met, all the components involved are previously sanitized and the procedure itself is validated according to ISO 17665-1 standard.

Unit: Steam sterilizer with fractionated vacuum and overpressure according to EN 13060 standard.

Procedure: Package each component to be treated in sterile barrier system or packaging complying with EN 11607 standard. Place the packed components in the steam sterilizer. Run the sterilization cycle according to the instructions manual of the device, selecting a temperature of 134° C and a duration of 10 minutes Storage: store the sterilized components as per instructions of the sterile barrier package or system chosen.

AIR FILTER

The device has a suction filter (A3) that should be replaced when it is dirty or changes color. Do not wash or reuse the same filter. Regular replacement of the filter is required to help and ensure correct performance of the compressor. Check the filter on a regular basis. Contact your local dealer or authorized service center for spare filters.

To replace the filter: Insert a flathead screwdriver between the edge of the filter and the body. Lift the filter and remove it by turning and pulling it upwards. The filter has been designed so that it is always fixed in its seat. **Use original FLAEM accessories only.**

SPARE PARTS

Description	Code
- RF6 Plus nebuliser kit	ACO234P
- 2 m connection tube	ACO225P
- Adult mask	ACO79
- Pediatric mask	ACO89
- Spare air filter kit (2 pcs)	ACO164P

SYMBOLS

TÜV Certification ref. EN 60601-1 3rd Ed. - EN 60601-1-11

 0051 CE Medical Marking ref. Dir. 93/42 EEC and subsequent amendments

Class II device

Important: check the operating instructions

“OFF” for part of equipment

“ON” for part of equipment

Complies with: European standard EN 10993-1 “Biological Evaluation of Medical Devices” and European Directive 93/42/EEC “Medical Devices”. Phthalates free. In conformity with Reg. (EC) no. 1907/2006

IP21

Enclosure protection rating: IP21.
(Protected against solid bodies over 12 mm. Protected against access with a finger;
Protected against vertically falling water drops.)

Manufacturer

Type BF applied part

Alternating current

Hazard: electrocution.
Consequence: Death.
Do not use device while taking a bath or shower

SN

Serial number of device

TECHNICAL SPECIFICATIONS

Mod. NEBULAIR

Voltage:	230V~ 50Hz 210VA	115V~ 60Hz 230VA	220V~ 60Hz 230VA	100V~ 50/60Hz
Approvals:				

Max Pressure:	3.5 ± 0.5 bar
Compressor air output:	14 LPM approx
Noise level (at 1 m):	55 dB (A) approx
Operation:	Continuous use
Dimensions	20(W)x30(D)x10(H) cm
Weight	2.400 Kg
Operating conditions:	Temperature: min 10°C; max 40°C Air humidity: min 10%; max 95%
Storage conditions:	Temperature: min -25°C; max 70°C Air humidity: min 10%; max 95%
Operating/storage atmospheric pressure:	min. 69KPa; máx. 106KPa

APPLIED PARTS

Type BF applied parts include: (C1,C2,C3,C4)

RF6 Plus NEBULIZER

Medication minimum capacity : 2 ml

Medication maximum capacity: 8 ml

Operating pressure bar	⁽²⁾ MMAD μm	⁽²⁾ BREATHABLE FRACTION < 5 μm (FPF) %	⁽¹⁾ DELIVERY ml/min	
			open valve-holder MAX	closed valve-holder MIN
1.2	2.8	76	0.54 approx	0.21 approx

(1) Data shown is in accordance with Flaem Nuova internal procedure I29-P07.5.

(2) In vitro testing certified by TÜV Rheinland Product Safety GmbH - Germany in compliance with the European Standard EN 13544-1. More details are available on request.

Disposal: In conformity with Directive 2012/19/EC, the symbol shown on the appliance to be disposed of indicates that this is considered as waste and must therefore undergo "sorted waste collection". The user must therefore take (or have taken) the above waste to a pre-sorted waste collection centre set up by the local authorities, or else give it back to the dealer when purchasing a new appliance of the same type. Pre-sorted waste collection and the subsequent treatment, recovery and disposal operations favour the production of appliances made of recycled materials and limit the negative effects of any incorrect waste management on the environment and public health. The unlawful disposal of the product by the user could result in administrative fines as provided by the laws transposing Directive 2012/19/EC of the European member state in which the product is disposed of.

ELECTROMAGNETIC COMPATIBILITY

This device has been designed to satisfy requirements currently required for electromagnetic compatibility (EN 60601-1-2:2007). Electrical medical devices require special care. During installation and use with respect to EMC requirements, it therefore required that they be installed and/or used according to the manufacturer's specification. Potential risk of electromagnetic interference with other devices, in particular with other devices for analysis and treatment. Radio and mobile telecommunications devices or portable RF (mobile phones or wireless connections) may interfere with the operation of electrical medical devices. For further information visit www.flaemnuova.it. Flaem reserves the right to make technical and functional changes to the product without notice.

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© 2015 FLAEM NUOVA®
FLAEM NUOVA S.p.A.