## Instructions for Use

### Cutimed® Siltec® B

#### Product description

Cutimed® Siltec® B dressings are an assortment of silicone coated foam dressings. Cutimed® Siltec® B is a sterile, single use, absorbent polyurethane foam dressing for atraumatic dressing changes that contains super-absorbent stripes which absorb and lock wound exudate. The wound contact surface is a perforated silicone layer that allows adherence to the periwound skin but not to the moist wound bed or to newly formed epithelial tissue. This minimises trauma and pain during dressing changes. The outer film is water-repellent, yet permeable to oxygen and vapor. The foam dressing and the super-absorbent stripes are designed to absorb and to lock away excess wound fluid inside the dressing, promoting a moist wound environment and minimizing the risk of maceration. Cutimed® Siltec® B is highly conformable and provides protection of the wound site. Cutimed® Siltec® B leaves no residues in the wound. The dressing reliably retains wound exudate, even under compression. If necessary the dressing can be cut to size. Cutimed® Siltec® B has an additional silicone adhesive border, that allows secure and gentle fixation of the dressing.

#### Intended purpose

Cutimed® Siltec® B is intended for the treatment of exuding wounds with low to high exudate levels such as venous and arterial ulcers, pressure ulcers, diabetic foot ulcers, surgical incisions, skin grafts and donor sites, lacerations or abrasions. Cutimed® Siltec® B may assist in the prevention of pressure ulcer as part of a comprehensive plan of continuous care, risk assessment and preventive care by a healthcare facility and healthcare professional.

Please contact your healthcare professional if you are unsure whether the product is appropriate for you.

#### Contraindications

There are no known contraindications.

#### Warnings and precautions

Cutimed® Siltec® B is packaged for single use. Do not re-use or re-sterilize as there is a risk of transmission of body fluids or contaminated tissue between patients. Do not use if the pouch is already open or damaged as the sterility of the device is guaranteed only when the pouch is unopened and undamaged prior to use. Discard open or unused material. Do not use with oxidizing solutions such as hypochlorite or hydrogen peroxide. For external use only. The wound should be inspected for signs of infection and treated according to clinical practice if required. In rare cases skin reactions (e.g. redness, itching) may occur.

If the treated condition deteriorates, fails to improve or if a side effect is observed, consult a physician or an appropriate healthcare professional.

#### Directions for use

- Clean the wound and your hands according to good clinical practice. If you are unsure how to clean the wound appropriately, please contact your healthcare professional for further instructions.
- Make sure, the surrounding skin is dry.
- Choose an appropriate dressing size that covers the wound and overlaps the wound edges. **(A)**
- Please ensure, that the pouch is not open or damaged. Open the pouch and take out the dressing.
- Remove protective film **(B)** and apply with the adherent side towards the wound. **(C)**
- Avoid stretching the dressing.
- Depending on wound location it may be appropriate to apply additional dressing fixation, e.g. in order to increase resistance to incontinence.
- Cutimed® Siltec® B may be used on an infected wound in addition to antimicrobial therapy under medical supervision.
- If necessary, Cutimed® Siltec® B can be used in combination with other dressings e.g. Cutimed® Sorbact®/Sorbact® (in case of infected wounds) or hydrogel e.g. Cutimed® Gel (for dry, necrotic wounds).
- Cutimed® Siltec® B can remain in place for up to 7 days, depending on wound conditions and saturation of wound dressing. It is recommended to change the dressing every 24 hours initially, moving to less frequent changes as necessary.
- Carefully remove any additional dressing fixation and gently remove Cutimed® Siltec® B from the wound and discard.
- If required, cleanse the wound prior to application of a new dressing. The dressing can be used in conjunction with radiotherapy. The dressing is showerproof.

#### Storage and Disposal

Store dressing away from direct sunlight at ambient temperature and humidity. Keep out of the reach of children. Discard used product according to local regulations, e.g. via regular household or clinical waste.

#### Notice

Any serious incident that has occurred in relation to this device should be reported to BSN medical GmbH and the competent authority of your state.

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Consult instructions for use  
Consulter les instructions d'usage

Medical device  
Dispositif médical

|                                                                                      |                            |                                            |                                                                                                                                |
|--------------------------------------------------------------------------------------|----------------------------|--------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------|
|    | Artwork ID:<br>JBR16984.01 | Review Loop:<br>10                         | <div>Printing Colours:</div> <div><div></div> PANTONE 287 C</div> <div><div></div> Magenta</div> <div><div></div> Yellow</div> |
| Project Title:<br>C2365_073_Cutimed_Siltec_B                                         |                            | Date of Review Loop:<br>30.11.2021 JR_KLND |                                                                                                                                |
| Contact<br>ArtworkManagementMedical@essity.com                                       |                            | Replaces Artwork ID:<br>JBR15313.04        |                                                                                                                                |
| Contact Address: BSN medical GmbH<br>Schützenstraße 1 – 3<br>22761 Hamburg   Germany |                            | Technical Drawing for print ID             |                                                                                                                                |
| Languages:<br>EN, FR, DE, NL, IT, ES, PT, SV, FI, NO, DA                             |                            | Note                                       |                                                                                                                                |

The following artwork was approved:

|                     |                                     |
|---------------------|-------------------------------------|
| Asset Name          | JBR16984.01                         |
| Component Name      | JBR16984.01                         |
| Project Title       | C2365_073 CM Siltec B / Sacrum IFUs |
| Revision            | 1                                   |
| Reason for Revision | new CM design                       |

## Asset Status

| Date (GMT)                 | Signed By           | Role                   | Action                                                               | Comment                                                 |
|----------------------------|---------------------|------------------------|----------------------------------------------------------------------|---------------------------------------------------------|
| 2021-11-18<br>11:41:47.0   | Ann-Kristin Wächter | Artwork Management     | Artwork Submitted for Approval                                       | I am the author of this document                        |
| 2021-12-01<br>08:50:08.052 | Kim Roettger        | Regulatory Affairs HAM | Artwork Approved : Electronic Signature (HMS Medical GM&I PM GRA MA) | I attest to the accuracy and integrity of this document |
| 2021-12-01<br>08:22:17.736 | Yvonne Link         | Product Development    | Artwork Approved : Electronic Signature (HMS Medical GM&I PM GRA MA) | I attest to the accuracy and integrity of this document |
| 2021-12-01<br>07:27:57.799 | Helgo Schieber      | Packaging Development  | Artwork Approved : Electronic Signature (HMS Medical GM&I PM GRA MA) | I attest to the accuracy and integrity of this document |
| 2021-12-01<br>07:22:51.938 | Stephanie Krause    | Scientific Affairs     | Artwork Approved : Electronic Signature (HMS Medical GM&I PM GRA MA) | I attest to the accuracy and integrity of this document |
| 2021-12-01<br>09:54:22.0   | Ann-Kristin Wächter | Artwork Management     | Final Sign Off                                                       | I attest to the accuracy and integrity of this document |

THIS REPORT IS CREATED FROM ELECTRONIC SIGNATURE ACTIVITY WITH KALLIK AMS 360

CREATED ON 2021-12-01 09:54:22.0