SINCE 1966
MADE IN ITALY

Plus RF6

by **FLAEM**®

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AEROSOL THERAPY NEBULIZER

Congratulations on your purchase and thank you for choosing our unit.

Our goal is to fully satisfy customers by offering them cutting-edge systems for the treatment of respiratory tract ailments.

Carefully read this instruction manual and store it for future reference. Only use the accessories as described in this manual. This is a medical device intended to nebulize and administer medications prescribed or recommended by a doctor.

Visit our Internet site www.flaem.it to view the whole range of Flaem products.

ACCESSORIES THAT CAN BE USED WITH THIS NEBULIZER

- 1 - Child nasal inhaler
- 2 - Adult nasal inhaler
- 3 - Child bi-material face mask
- 4 - Adult bi-material face mask
- 5 - Polyethylene (PE) face mask
- 6 - Medical PVC child face mask
- 7 - Medical PVC adult face mask
- 8 - Mouthpiece with valve
 - 8.1 - Expiratory valve
- 9 - Non-invasive nasal prong
- 10 - Baby mask set (see its instruction manual)
- 11 - Nebulisation manual control
- 12 - Connector
- 13 - 1 or 2 m Connection hose
- 14 - RF6 Plus with ambient contamination prevention set (see its instruction manual)

NEBULIZER

- 15 - RF6 Plus NEBULISER
 - 15.1 - Lower part
 - 15.2 - Complete nozzle
 - 15.3 - Upper part with valve-holder
 - 15.4 - Double action inspiratory valve

WARNINGS

- This is a medical device (compliant with Dir. 93/42/EEC) and may only be used with medication prescribed or recommended by a doctor. Before using the accessories, consult the instruction manual of the device. It is important that the patient reads and understands the information on how to use the device and the nebulizer. Contact your retailer or nearest service center for any questions.
- When used very frequently, we recommend that the nebuliser be replaced every 6 months (or earlier if it becomes clogged), in order to guarantee the maximum therapeutic effect.
- Children and disabled persons should always use the nebulizer under strict supervision of an adult who has read this manual.
- The nebulizer contains small parts that may be swallowed by children; keep therefore the nebulizer out of the reach of children

OPERATING INSTRUCTIONS

Before each use, wash your hands accurately and clean the device as described in the paragraph "CLEANING SANITISATION DISINFECTION STERILISATION".

The nebuliser cup and accessories are for personal use only to avoid potential risks of transmitting infection.

During application, it is advisable to suitably protect yourself from any dripping.

1. Assemble the nozzle as indicated in the "Connection diagram" in point 15.2 Insert the complete nozzle in the bottom part (15.1). Insert the Double action inspiratory valve (15.4) in the upper part (15.3) as shown in the "Connection diagram" in point 15.
2. Pour the medication prescribed by the physician into the lower part (15.1). Close the

- nebuliser by turning the upper part (15.3) clockwise.
3. Connect the accessories as indicated in the "Assembly diagram." Sit comfortably holding the nebuliser in your hand, place the mouthpiece onto your mouth or alternatively use the nosepiece, mask or newborn set.
 4. Start the device (see its operating instructions). Gently take a deep breath; we recommend holding your breath for an instant to allow the inhaled drops of aerosol to deposit. Then exhale slowly.
 5. Once the process is complete, switch off the device

ATTENTION: If after the therapy session an evident deposit of moisture forms within the pipe (13), detach the pipe from the nebuliser and dry it with ventilation from the said compressor; this operation prevents possible proliferation of mould inside the pipe.

USAGE INSTRUCTIONS FOR THE "RF6 PLUS" NEBULISER WITH VALVE SYSTEM AND DUAL VELOCITY FUNCTION.

Thanks to the design of the RF6 Plus interior flask ducts, it has been possible to obtain the right granulometry, which allows for the treatment to reach even the lower respiratory tract. RF6 Plus is indicated for the administration of all kinds of medications, even the most costly, and even to patients with chronic diseases. The effective valve system, combined with the medication anti-dispersal saving device, allows for high therapeutic efficacy by maximizing the intake of the medication in the lungs and minimizing its dispersal in the room. The RF6 Plus is particularly indicated for pediatric use, as the non-spill system prevents the medicine from being spilt, even by accidental overturning (its efficacy has been tested up until 5 ml, the maximum medication dose in the most common therapies). The nebuliser makes it possible to deliver the medication in 2 ways:


CLOSED VALVE-HOLDER

In this case the inhalation therapy allows the maximum drug intake and therefore the maximum therapeutic effect, minimizing the loss of medication. This configuration is to be used with acces-


sories 8 - 10 - 14.

OPEN VALVE-HOLDER


For a quicker treatment. This configuration is to be used with accessories 1 - 2 - 3 - 4 - 5 - 6 - 7 - (8 with nasal inhaler 9).

SOFTTOUCH FACE MASKS


Soft biocompatible material
Dispersion Limiting Vent

SoftTouch masks have an outer edge made of **soft biocompatible material** that ensures excellent adherence to the face, and is also equipped with an **innovative Dispersion Limiting Device**. These distinctive elements that distinguish it allow greater sedimentation of medication in the patient and also **limit dispersion**.


During the inspiratory phase, the **tab**, acting as a **Dispersion Limiting Device**, bends inwards towards the mask.


During the expiratory phase, the **tab**, acting as a **Dispersion Limiting Device**, bends outwards from the mask.

USE OF NEBULIZATION MANUAL CONTROL

To achieve continuous nebulisation action you should not use the manual nebulisation control (11), especially in the case of children or persons with reduced physical, sensory, or mental capabilities. The manual nebulisation control is useful for limiting dispersion of the medication in the surrounding environment.

To start nebulizing close with a finger the hole of the nebulizer manual control (11) and breathe in gently; we recommend to hold your breath for a moment so that the inhaled aerosol droplets can be depos-

meanwhile, to disable nebulizing, remove your finger from the hole of the nebulizer manual control to avoid the waste of drug, optimizing its acquisition. Then exhale slowly.

ited,

CLEANING SANITISATION DISINFECTION STERILISATION

Switch off the device before any cleaning procedure and unplug the power cable from the socket.

DEVICE AND TUBING EXTERIOR

Use only a damp cloth with antibacterial soap (non-abrasive and with no solvents of any sort).

ACCESSORIES

Open the nebuliser by turning the upper part (15.3) anticlockwise, remove the nozzle (15.2) from the upper part (15.3). Disassembly by pressing as shown by the 2 arrows in the "Connection diagram" in point 15.2.

Then proceed according to the following instructions.

CLEANING AT HOME - SANITISATION AND DISINFECTION

SANITISATION

Before and after each use, sanitise the nebuliser cup and the accessories, choosing one of the methods described below.

method A: Sanitise accessories **1-2-3-4-5-6-7-8-9-11-12-15.1-15.2-15.3-15.4** under potable hot water (approximately 40°C) with a gentle, non abrasive dish detergent.

method B: Sanitise accessories **1-2-3-4-5-6-7-8-9-11-12-15.1-15.2-15.3-15.4** in the dishwasher with a hot cycle.

method C: Sanitise accessories **1-2-3-4-5-6-7-8-9-11-12-15.1-15.2-15.3-15.4** by immersing them in a solution of 50% water and 50% white vinegar, then rinse thoroughly under potable hot water (approximately 40 °C).

If you want to also perform the cleaning for DISINFECTION, jump to the DISINFECTION paragraph.

After having sanitised the accessories, shake them vigorously and place them on a paper towel. Alternatively, dry them with a jet of hot air (for example, a hair dryer).

DISINFECTION

After sanitising the nebuliser cup and the accessories, disinfect them choosing one of the methods described below.

method A: Accessories **1-2-3-4-5-6-7-8-9-11-12-15.1-15.2-15.3-15.4** can be disinfected. The disinfectant must be an electrolytic chloroxidizer (active principle: sodium hypochlorite) specific for disinfecting, which is available in any pharmacy.

Implementation:

- Fill a container big enough to hold all of the parts to disinfect with a solution of potable water and disinfectant, according to the proportions indicated on the packaging of the disinfectant.
- Completely immerse each part in the solution, taking care to avoid the formation of air

- bubbles on the parts. Leave the parts immersed for the amount of time indicated on the packaging of the disinfectant associated with the concentration chosen for the solution.
- Remove the disinfected parts and rinse abundantly with warm potable water.
 - Dispose of the solution following the instructions provided by the disinfectant manufacturer.

method B: Disinfect the accessories **1-2-3-4-8-9-11-12-15.1-15.2-15.3-15.4** by boiling them in water for 10 minutes; use demineralised or distilled water to prevent calcium deposits.

method C: Disinfect the accessories **1-2-3-4-8-9-11-12-15.1-15.2-15.3-15.4** with a hot steam steriliser for baby-bottle (not the microwave type). Perform the process faithfully following the instructions of the steriliser. To ensure that the disinfection is effective, choose a steriliser with an operating cycle of at least 6 minutes.

After having disinfected the accessories, shake them vigorously and place them on a paper towel. Alternatively, dry them with a jet of hot air (for example, from a hair dryer).

CLEANING IN A CLINICAL OR HOSPITAL SETTING - DISINFECTION AND STERILISATION

Before disinfection or sterilisation, sanitise the nebuliser cup and the accessories, choosing one of the methods described below.

method A: Sanitise accessories **1-2-3-4-5-6-7-8-9-11-12-15.1-15.2-15.3-15.4** under potable hot water (approximately 40°C) with a gentle, non abrasive dish detergent.

method B: Sanitise accessories **1-2-3-4-5-6-7-8-9-11-12-15.1-15.2-15.3-15.4** in the dishwasher with a hot cycle.

DISINFECTION

Accessories **1-2-3-4-5-6-7-8-9-11-12-15.1-15.2-15.3-15.4** can be disinfected.

The disinfectant must be an electrolytic chloroxidizer (active principle: sodium hypochlorite) specific for disinfecting, which is available in any pharmacy.

Implementation:

- Fill a container big enough to hold all of the parts to disinfect with a solution of potable water and disinfectant, according to the proportions indicated on the packaging of the disinfectant.
- Completely immerse each part in the solution, taking care to avoid the formation of air bubbles on the parts. Leave the parts immersed for the amount of time indicated on the packaging of the disinfectant associated with the concentration chosen for the solution.
- Remove the disinfected parts and rinse abundantly with warm potable water.
- Dispose of the solution following the instructions provided by the disinfectant manufacturer.

If you want to also perform the STERILISATION, jump to the STERILISATION paragraph.

After having disinfected the accessories, shake them vigorously and place them on a paper towel. Alternatively, dry them with a jet of hot air (for example, from a hair dryer).

STERILISATION

Accessories **1-2-3-4-8-9-12-15.1-15.2-15.3-15.4** can be sterilised.

Device: Fractionated vacuum overpressure steam steriliser in accordance with EN 13060.

Implementation: Wrap every single part to be treated with a sterile barrier system or packaging in accordance with Norm EN 11607. Place the packed components in the steam steriliser. Run the sterilisation cycle according to the operating instructions of the device by selecting a temperature of 134 °C and a time of 10 minutes first.

Storage: Store the sterilised parts as per the instructions for use of either the sterile barrier system or packaging.

The sterilisation procedure is validated in its conformity to ISO 17665-1.

For the “**CLEANING SANITISATION DISINFECTION STERILISATION**” of accessories 10-14

follow the instructions in the instruction manual.

TECHNICAL SPECIFICATIONS

RF6 Plus Nebulizer

In compliance with 93/42/EEC Directive

Minimum fill volume: 2 ml

Maximum fill volume: 8 ml

Operating conditions:

Temperature: min 10°C; max 40°C

Air humidity: min 10%; max 95%

Atmospheric pressure: min 69KPa; max 106KPa

Storage conditions:

Temperature: min -25°C; max 70°C

Air humidity: min 10%; max 95%

Atmospheric pressure: min 69KPa; max 106KPa

APPLIED PARTS

The applied parts of type BF are:

patient accessories (1, 2, 3, 4, 5, 6, 7, 8, 9, 10, 14)

FLAEM compressor range	Operating pressure	MMAD (μm) ⁽²⁾	BREATHABLE FRACTION < 5 μm (FPF) % ⁽²⁾	DELIVERY ml/min ⁽¹⁾ (approx)	
				open valve-holder MAX	closed valve-holder MIN
F700 / F1000	0.8 bar	3.2	73	0.42	0.26
F2000	1.2 bar	2.8	76	0.54	0.31

⁽¹⁾ Data shown is in accordance with Flaem Nuova internal procedure I29-P07.5. ⁽²⁾ In vitro characterization carried out at the Inamed Research GmbH Co. KG on behalf of TÜV Rheinland LGA Products GmbH – Germany in compliance with the new European Standard for aerosol therapy units. Standard EN 13544-1. Further details are available upon request.

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