Super Absorbent Hydrogel Colloidal Sheet

Hoch absorbierende kolloidale Hydrogeleplatte

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PRODUCT DESCRIPTION

XTRASORB[®] HCS Wound Dressing is a hybrid dressing combining hydrogel and hydrocolloid technology, with an additional super absorbent polymer for added fluid handling capability. The dressing consists of a water/glycerin component to help hydrate dry to lightly exuding wounds, and a sodium carboxymethylcellulose (CMC) /super absorbent polymer (SAP) component to help absorb moderate amounts of exudates. On dry to lightly exuding wounds, the hydrating dressing is cooling and soothing upon application. On moderately to exuding wounds, the CMC/SAP forms a hydrophilic gel that locks exudates into the dressing.

XTRASORB[®] HCS is available in the form of both a borderless and a bordered "island" dressing, both backed by an integral PU film barrier that provides a barrier to external contaminants. On the bordered version, this PU film extends beyond the center wound pad, creating a skin-friendly adhesive border to aid fixation. The dressing has been formulated and designed to encourage wound bed preparation, granulation and subsequent epithelialization of chronic wounds, while minimizing the risk of infection by acting as a bacterial barrier. The dressings help to provide an environment for optimally balanced moist wound healing by protecting against wound dehydration while also absorbing excess wound exudates. The unique gelling action of the CMC/SAP is designed to lock wound fluid within the dressing, helping to reduce the risk of maceration.

INDICATIONS

Under the supervision of a health care professional, XTRASORB[®] HCS dressing may be used for the management of acute and chronic, non-infected or mildly infected wounds. It may be used on dry, lightly exuding, and moderately exuding wounds including:

- Pressure Ulcers
- Venous leg Ulcers
- Arterial Ulcers
- Diabetic Foot Ulcers
- Post-Op Wounds
- Traumatic Wounds
- Donor Sites

CONTRAINDICATIONS

Known hypersensitivity to the product itself or to its components, including glycerin and sodium carboxymethylcellulose.

Not for use on heavily bleeding wounds or third degree burns.

PRECAUTIONS AND OBSERVATIONS

Sterility of the dressing is guaranteed unless the pouch is damaged or opened prior to use. DO NOT use if either is occurred.

The product is intended for single use only.

Do not use after expiry date.

In the case of intolerance to the dressing, remove and clean the area carefully. During the body's natural healing process, non-viable tissue is removed from the wound (autolytic debridement), which could initially make the wound appear larger. If the wound continues to grow larger after the first few dressing changes, consult a healthcare professional. The wound should be inspected during dressing changes. Consult a healthcare professional if there are signs of infection (increased pain, increased redness, wound drainage); excessive bleeding, an unexpected change in wound color and/or odor; irritation (increased redness and/or inflammation); sensitivity (allergic reaction); no signs of healing.

If difficulty is experienced on removing the dressing, it should be irrigated or soaked with water or a sterile saline solution.

Because this dressing provides an environment that supports the growth of new blood vessels, the delicate newly-formed tissue may produce blood stained fluid.

Appropriate supportive measures should be taken where necessary (e.g. use of graduated compression in the management of venous leg ulcers or pressure relief measures in the management of pressure ulcers, systemic antibiotics and frequent monitoring in the treatment of wound infection, control of blood glucose for diabetic ulcers, etc.)

A change in dressing regimen can result in an initial increased level of exudate, which temporarily may require an increased dressing change frequency.

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