Gebrauchsinformation

DE

KD intraartikulär Gel

Natriumhyaluronat intraartikulär

Formel	KD intraartikulär Gel					
1 ml enthält:	1.0%	1.2%	1.4%	1.6%	2.2%	3.0%
Natriumhyaluronat	10.000 mg	12.000 mg	14.000 mg	16.000 mg	22.000 mg	30.000 mg
NaCl 0.9%	8.500 mg	8.500 mg	8.500 mg	8.500 mg	8.500 mg	8.500 mg
Dinatriumhydrogenphosphat	0.563 mg	0.563 mg	0.563 mg	0.563 mg	0.563 mg	0.563 mg
Natriumhydrogenphosphat	0.045 mg	0.045 mg	0.045 mg	0.045 mg	0.045 mg	0.045 mg
Wasser für Injektionszwecke	q.s.	q.s.	q.s.	q.s.	q.s.	q.s.

Quelle

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Origin

Sodium hyaluronate is obtained by bacterial fermentation and is used as a viscoelastic gel to supplement the intra-articular connective tissue. It is sterile and apyrogen.

Composition

1 ml of solution contains 10 mg, 12 mg, 14 mg, 16 mg, 22 mg or 30 mg of sodium hyaluronate with fermentative origin, a low phosphate tampon sterile, isotonic and apyrogen, pH 7.2 and water for injection.

Indications

It is injected into the knee to relieve pain and decreased articular mobility associated with degenerative knee lesions and other synovial joints.

Contra-Indications

The product should not be used for patients who developed a known hypersensitivity reaction to any of the components of the product.

Cautions

The product should be injected into the articulation. Intra-vascular or intra-tissular injection should be avoided.

Since no clinical data are available on the use of hyaluronic acid in children or pregnant woman, these patients should not be treated.

If packaging is damaged product should not be used.

The product has to be used before the expiry date indicated on the packaging.

No sign of bioincompatibility has been detected throughout biocompatibility testing of the product. Some local transient reactions may occur after intra-articular injection, such as redness, swelling, heat and pain. The local application of ice for 5-10 minutes may reduce these effects.

No known incompatibility with other intra-articular injectable products has been reported up to now.

Simultaneous intake of analgesics or anti inflammatory drugs by oral route may be useful during the first days after treatment.

The product must not be used in association with quarternary ammonium compounds.

The product must not be resterilised or reuse because this could lead to a risk of infection and/or cross-infection and/or malfunction of the product.

Administration and dosage

The syringe is prefilled, the Luer lock closure is removed and a suitable sharp edged tip sterile needle is screwed on the Luer lock. Use of a single use needle is recommended. On desire, needles are included. Remove air bubbles from the syringe before injection.

The skin is disinfected with local antiseptic.

The product is injected in the joint, once a week, 3-5 times. It is recommended to renew the injection, about 2 weeks later. The beneficial effects of 3-5 sequential treatments are maintained during at least 6 months. In the case of intra articular swelling, with fluid should be removed by aspiration before the hyaluronic acid injection.

Properties and mode of action

The synovial fluid contains physiologically hyaluronic acid which is responsible for viscoelastic properties allowing, through lubricant properties, painless joint movements.

It also participates to the nutrition of the cartilage.

Intra-articular supplementation of the synovial fluid by hyaluronic acid injections improves the viscoelastic properties of the synovial fluid. After 3-5 intra-articular injections, with joint mobility is improved, and the pain associated with degenerative lesions is relieved.

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ALBOMED GmbH