Pulsoximeter PO-300 - Bedienungsanleitung
Pulse Oximeter PO-300 - Instruction manual
Pulsioxímetro PO-300 - Manual de instrucciones
Oxymètre de Pouls PO-300 - Mode d'emploi
Pulsossimetro PO-300 - Istruzioni per l'uso
Pulsoximeter PO-300 - Gebruiksaanwijzing
Pulsoksymetr PO-300 - Instrukcja obsługi
Pulsoximeter PO-300 - Bruksanvisning

www.pulox.de

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Instructions to User

Dear users, thank you very much for purchasing the Pulse Oximeter.

This Manual is written and compiled in accordance with the council directive MDD93/42/EEC for medical devices and harmonized standards. In case of modifications and software upgrades, the information contained in this document is subject to change without notice.

The Manual describes, in accordance with the Pulse Oximeter's features and requirements, main structure, functions, specifications, correct methods for transportation, installation, usage, operation, repair, maintenance and storage, etc. as well as the safety procedures to protect both the user and equipment. Refer to the respective chapters for details.

Please read the User Manual carefully before using this product. The User Manual which describes the operating procedures should be followed strictly. Failure to follow the User Manual may cause measuring abnormality, equipment damage and human injury. The manufacturer is NOT responsible for the safety, reliability and performance issues and any monitoring abnormality, human injury and equipment damage due to users' negligence of the operation instructions. The manufacturer's warranty service does not cover such faults. Owing to the forthcoming renovation, the specific products you received may not be totally in accordance with the description of this User Manual. We apologize for this.

This product is a medical device, which can be used repeatedly.

Warning

- ! Uncomfortable or painful feeling may appear if using the device ceaselessly, especially for the microcirculation barrier patients. It is recommended that the sensor should not be applied to the same finger for over 2 hours.
- ! For special patients, there should be a more prudent inspection in the placing process. The device can not be clipped on the edema and tender tissue.
- ! The light (infrared is invisible) emitted from the device is harmful to the eyes, so the user and maintenance personnel should not stare at the light.
- ! Testee can not use enamel or other makeup.
- ! Testee's fingernail can not be too long.
- ! Please refer to the correlative literature about the clinical restrictions and caution.
- ! This device is not intended for diagnosis.

The User Manual is published by our company. All rights reserved.

1. Safety

1.1 Instructions for Safe Operations

- ◇ Check the main unit and all accessories periodically to make sure that there is no visible damage that may affect patient's safety and monitoring performance about cables and transducers. It is recommended that the device should be inspected once a week at least. In case of obvious damage, stop using the device.
- ◇ Necessary maintenance must be performed by qualified service engineers ONLY. Users are not permitted to maintain it by themselves.
- ◇ The oximeter cannot be used together with devices not specified in User's Manual. Only accessory appointed or recommended by manufacture can be used with this device.
- ◇ This product is calibrated before leaving factory.

1.2 Warning

- ! Explosive hazard—DO NOT use the oximeter in environment with inflammable gas such as some ignitable anesthetic agents.
- ! DO NOT use the oximeter while the testee measured is being measured by MRI and CT.
- ! Be careful with the use of the lanyard cord. Improper use of the lanyard cord will cause device damage not covered under the manufacturer's warranty. Swinging the device by the lanyard cord will void the warranty. Please do not use lanyard cord if allergic to lanyard cord.
- ! The disposal of scrap instrument and its accessories and packings (including battery, plastic bags, foams and paper boxes) should follow the local laws and regulations.
- ! Please check the packing before use to make sure the device and accessories are totally in accordance with the packing list, or else the device may have the possibility of working abnormally.
- ! Please choose the accessories and probe which are approved or manufactured by the manufacturer, or else it may damage the device.
- ! Please choose the battery chargers which should be ensured compliance with the requirements of IEC 60601-1, or else it may damage the device.
- ! Please don't use the device in the course of charging.
- ! Please don't measure this device with functional tester for the device's related information.

1.3 Attention

- Keep the oximeter away from dust, vibration, corrosive substances, explosive materials, high temperature and moisture.
- If the oximeter gets wet, please stop operating it.
- When it is carried from cold environment to warm or humid environment, please do not use it immediately.
- DO NOT operate keys on front panel with sharp materials.
- High temperature or high pressure steam disinfection of the oximeter is not

permitted. Refer to User Manual in the relative chapter (7.1) for instructions of cleaning and disinfection.

- Do not immerse the oximeter in liquid. When it needs cleaning, please wipe its surface with disinfectant solution using a soft cloth. Do not spray any liquid directly onto the device.
- When cleaning the device with water, the temperature should be lower than 60°C.
- As to the fingers which are too thin or too cold, it would probably affect the normal measure of the patients' SpO₂ and pulse rate, please clip a thicker finger such as thumb and middle finger deeply enough into the probe.
- The update period of data is less than 5 seconds, which is changeable according to different individual pulse rate.
- Please read the measured value when the waveform on screen is equally and steady-going, This measured value is optimal value. And the waveform at the moment is the standard one.
- If some abnormal conditions appear on the screen during test process, pull out the finger and reinsert to restore normal use.
- The device has an estimated lifespan of three years from the first electrified use.
- The hanging rope attached to the device is made from Non- allergy material, if particular group are sensitive to the hanging rope, stop using it. In addition, pay attention to the use of the hanging rope, do not wear it around the neck on the purpose of avoiding harm to the patient.
- This device includes an alarmfunction, users can check this according to chapter 6 as a reference.
- The device has the function of limits alarming, when the measured data is beyond the highest or lowest limit, the device would start alarming automatically on the premise of the alarming function is on.
- The device may not work for all patients. If you are unable to achieve stable readings, discontinue use.
- A flexible circuit connects the two parts of the device. Do not twist or pull on the connection.

2. Overview

The pulse oxygen saturation is the percentage of HbO₂ in the total Hb in the blood, in other words the O₂ concentration in blood. It is an important bio-parameter for respiration. A number of diseases relating to the respiratory system may cause a decrease of SpO₂ in the blood, furthermore, some other causes such as the malfunction of human body's self-adjustment, damages during surgery, and the injuries caused by some medical checkup would also lead to the difficulty of oxygen supply in human body, and the corresponding symptoms would appear as a consequence, such as vertigo, impotence, vomit etc. Serious symptoms might bring danger to a person's life. Therefore, prompt information of patients' SpO₂ is of great help for the doctor to discover the

potential danger, and is of great importance in the clinical medical field. The Pulse Oximeter features in small volume, low power consumption, convenient operation and being portable. It is only necessary for patients to put one of his fingers into a probe for diagnosis, and a display screen will directly show the measured value of pulse oxygen saturation with the high veracity and repetition.

2.1 Features

- A. Operation of the product is simple and convenient.
- B. The product is small in volume, light in weight and convenient for carrying.
- C. Low power consumption

2.2 Major Applications and Scope of Application

The Pulse Oximeter can be used in measuring the pulse oxygen saturation and pulse rate through finger. The product is suitable for being used in family, hospital, oxygen bar, community healthcare, physical care in sports (It can be used before or after doing sports and it is not recommended to use the device during the process of exercising) and more.

 **The problem of overrating would emerge when the patient is suffering from toxicosis caused by carbon monoxide, the device is not recommended to be used under this circumstance.**

2.3 Environment Requirements

Storage Environment

- a. Temperature : $-40^{\circ}\text{C} \sim +60^{\circ}\text{C}$
- b. Relative humidity : $\leq 95\%$
- c. Atmospheric pressure : $500\text{hPa} \sim 1060\text{hPa}$

Operating Environment

- a. Temperature: $10^{\circ}\text{C} \sim 40^{\circ}\text{C}$
- b. Relative Humidity : $\leq 75\%$
- c. Atmospheric pressure: $700\text{hPa} \sim 1060\text{hPa}$

3. Principle

The principle of the Oximeter is as follows: An experience formula of data process is established taking use of Lambert Beer Law according to Spectrum Absorption Characteristics of Reductive Hemoglobin (Hb) and Oxyhemoglobin (HbO₂) in glow & near-infrared zones. The operation principle of the device is: Photoelectric Oxyhemoglobin Inspection Technology is adopted in accordance with Capacity Pulse Scanning & Recording Technology, so that two beams of different wavelengths of lights can be focused onto a human nail tip through perspective clamp finger-type sensor. The measured signal can be obtained by a photosensitive element, information acquired through which will be shown on screen through treatment in electronic circuits and microprocessor.


Figure 1

4. Technical Specifications

4.1 Main Performance

- A. Display of SpO₂ value, PR value, bar graph and pulse waveform.
- B. Low-battery indication: low-battery indication appears when the battery voltage is too low to work. Automatic standby function.
- C. Adjustable screen brightness, the display direction can be changed automatically.
- D. PR sound indication; sound prompt for over-limit, finger-out and low battery.
- E. Memory function
- F. The data can be uploaded to the terminal equipment by wired mode.
- G. Charging function
- H. Time display

4.2 Main Parameters

A. Measurement of SpO₂

Measuring range: 0%~100%

Accuracy: 70% ~ 100%, $\pm 2\%$; 0% ~ 69%, unspecified

B. Measurement of pulse rate

Measuring range: 30bpm~250bpm

Accuracy: ± 2 bpm or $\pm 2\%$ (select larger)

C. Resolution

SpO₂: 1%, Pulse rate: 1bpm

D. Measurement Performance in weak filling condition

SpO₂ and pulse rate can be shown correctly when pulse-filling ratio is 0.4%.

SpO₂ error is $\pm 4\%$, pulse rate error is ± 2 bpm or $\pm 2\%$ (select larger).

E. Resistance to surrounding light

The deviation between the value measured in the condition of man-made light or indoor natural light and that of darkroom is less than $\pm 1\%$.

F. Power supply requirement

3.6 V DC ~ 4.2V DC

G. Optical Sensor

Red light (wavelength is 660nm, 6.65mW)

Infrared (wavelength is 905nm, 6.75mW)

H. Adjustable alarm range

SpO₂: 0%~100%

Pulse Rate: 0bpm~254bpm

5. Installation

5.1 View of the front panel

Fig. 2 Appearance

USB interface: connect with USB cable

Button: exit standby, pause sound prompt, enter menu, menu operation.

RESET hole: reset button inside of it, be used to reset the device.

5.2 Interface introduction


Fig. 3 Measurement interface


Fig. 4 Sketch map for installing the hanging rope

Bluetooth sign (Bluetooth wired equipment) (see Figure 3).

5.3 Installing the hanging rope

A. Put the thinner side of the rope through the hole (see Figure 4).

B. Put the wider side of the rope through the thinner side which has been put

through the hole, then tighten it.

5.4 USB port

It is used to connect a personal computer to export the trend data or charge the lithium battery via a data line (see Figure 2).

5.5 Structure and accessories

- A. Structure: main unit, USB cable, power adapter
 B. Accessories: one hanging rope, one power adapter, one USB cable, PC software (downloadable), one User Manual

Please check the device and accessories according to the list to avoid that the device can not work normally.

6. Operating

6.1 Application method

6.1.1 Measurement

- A. Pinch the clamp, put a finger into the rubber hole, then release it, as shown in Figure 5.


Fig. 5 Sketch map for finger placement


Fig. 6 Time interface

- ⚠ When inserting the finger, the light emitting from the sensor must be directly irradiated to the side of the fingernail.
- ⚠ During measuring, do not shake the finger and keep quiet, not move.
- B. Press "Button" to exit from the standby mode, then it displays the measurement interface.
- ⚠ Please synchronize the time with the master device when using it for the first time, refer to chapter 6.1.3 for relative operations.

- C. Read the results directly from the screen after a period of time.
- D. Under measurement interface and non-alarm state, short press the button to enter the clock interface (as Figure 6).
- E. Under non-measurement state, it will enter standby mode automatically when there is no operation within 5s.

6.1.2 Pause sound prompt

- A. Sound prompt, including: over-limit, low-battery and finger out.
- B. Under the measurement interface, turn on the sound prompt, when the sound prompt occurs, short press the button to pause the sound prompt, and it will resume automatically after about 60s.
- C. If you want to turn off the sound prompt permanently, please set it in menu.

6.1.3 Menu operation

Under the measurement interface, long press the button to enter the main menu interface as shown in Figure 7, sound, record, clock, Bluetooth, etc. can be set, methods are as followings:


Fig. 7 Main menu


Fig. 8 Sound prompt menu


Fig. 9 Record menu

A. Sound prompt setting

Under main menu, short press the button to select "Sound", then long press the button to enter its setting interface shown in Figure 8.

Short press the button to select the option to be adjusted, then long press the button to change the value.

"Direction": direction, "up": increase the value, "down": decrease the value

"SpO₂ ALM HI": upper limit prompt for SpO₂ over-limit

"SpO₂ ALM LO": lower limit prompt for SpO₂ over-limit

"PR ALM HI": upper limit prompt for PR over-limit

"PR ALM LO": lower limit prompt for PR over-limit

"Alarm": prompt for over-limit, "off": close, "on": open.

"Pulse Sound": PR sound, "off": close, "on": open.

Lower limit can not exceed the upper limit, and the upper limit can not be lower than the lower limit when adjusting the values. SpO2 range: 0 % ~ 100 %, PR range: 0 ~ 254 bpm

The values displayed in Figure 8 are the initial values of over-limit prompt. After setting, short press the button to select "Exit", then long press the button to exit sound setting interface, and return to Main Menu interface.

B. Data storage

Under the main menu, short press the button to select "Record", then long press the button to enter the Record Menu interface as shown in Figure 7. It indicates that the device is storing when the red dot "R●" in measurement interface flickers.

Short press the button to select the option to be adjusted, then long press the button to change the value.

"Mode": record mode selection, including: "Auto" and "Manual" mode. Under "Manual" mode, select to turn on / off memory by "Record".

Auto record: start recording after stable data appear, pull out the finger to finish recording a group of data (99 group of data at most), the total duration does not exceed 72 hours.

Manual record: store up to 24-hour data.

When the memory is full, it will display "Memory is full!", then it will enter the standby mode after several seconds. When exiting the standby mode, it will display "Memory is full!" to give user a prompt, press the button again to enter the measurement interface.

Note: under manual mode, when "Record" is "ON", the device will prompt to clear the data stored last time.

It will display "Recording" when there is no operation under record state for 30s, then it will enter energy saving mode after several seconds, long press the button to exit this mode; short press the button, it will display "Recording".

Note: under data recording state, after the display screen turns off automatically, in order to save power, pulse sound indication will turn off automatically.

"Seg": data segment.

After setting, short press the button to select "Exit", long press the button to exit record menu and return to main menu.

"Delete All": delete all records (auto record mode is shown as Figure 9).

Note: please upload data in time after recording, otherwise the data may be covered when the storage space is full.

Note: the historical data will be deleted once switching the mode. Under recording state, the record mode can not be switched; under manual mode, the record mode can be switched only when turning off recording firstly.

C. Time setting

a. Connect the master device to synchronize device time

After connecting real-time transmission (refer to relative chapter (6.1.4) for the connection method), select "Options" - "Synchronize Device Time" on PC

software interface to synchronize the device time.

b. Set device time manually

Under main menu, short press the button to select “Clock”, long press the button to enter its sub-menu as shown in Figure 10.

Short press the button to select the option to be adjusted, then long press the button to change the value. The device adopts 24-hour clock.

“set time”: set the time, “yes”: to set the time, “no”: not set the time

“set year”: set the year

“set month”: set the month

“set day”: set the day

“set hour”: set the hour

“set minute”: set the minute

Adjustable range for year: 2015 ~ 2045, month: 1 ~ 12, day: 1 ~ 30 (when there are 31 days in a month, it is 1 ~ 31), hour: 1 ~ 23, minute: 1 ~ 59.

After setting, short press the button to select “Exit”, then long press the button to exit time setting interface and return to main menu.

Set Time	no
Set Year	2019
Set Month	01
Set Day	01
Set Hour	03
Set Minute	00
Exit	

Fig. 10 Time setting interface

System Menu	
Hard. Ver.	2.0.0
Soft. Ver.	2.0.2
ID	user
Demo	off
Sound Volume	3
Brightness	1
Exit	

Fig. 11 System menu interface

D. System setting and other options introduction

Under main menu, short press the button to select “System”, then long press the button to enter the interface as shown in Figure 11.

Short press the button to select the option to be adjusted, then long press the button to change the value.

“Hard.Ver.”: hardware version

“Soft.Ver.”: software version

“ID”: user name

“Demo”: set the Demo mode, “on”: turn on the Demo mode, “off”: turn off the Demo mode.

“Sound Volume”: set the sound volume, adjustable range: 1 ~ 3

“Brightness”: set the screen brightness, adjustable range: 1 ~ 4

After setting, short press the button to select “Exit”, then long press the button to exit the system menu and return to main menu.

E. Exit main menu

Under main menu, short press the button to select “Exit”, then long press the button to exit the main menu and return to the measurement interface.

6.1.4 Data upload

Connect the device to the computer by the USB cable, upload the data after connecting the PC software properly, refer to “Software operating instruction”

for details.

Note: the PC software can be downloaded from our official website.

6.1.5 Charging

Power adapter can be selected to charge for the device. It indicates that the device is charging when the indicator is light, the charging is finished when the indicator is off.

6.1.6 Reset

Pinch the clip and the finger pad in accordance with Figure 12, then use a pointed and hard object (for example, a paper clip) to press the reset button inside of the reset hole, to reset the device.


Figure 12 Sketch map for reset operation

6.2. Attention for operation

- A. Please check the device before using, and confirm that it can work normally.
- B. The finger should be in a proper position (see the attached illustration of Figure 5 for reference), or else it may result in inaccurate measure.
- C. The SpO₂ sensor and photoelectric receiving tube should be arranged in a way with the subject's arteriole in a position there between.
- D. The SpO₂ sensor should not be used at a location or limb tied with arterial canal or blood pressure cuff or receiving intravenous injection.
- E. Do not fix the SpO₂ sensor with adhesive or else it may result in venous pulsation and inaccurate measure of SpO₂ and pulse rate.
- F. Excessive ambient light may affect the measuring result. It includes fluorescent lamp, dual ruby light, infrared heater, direct sunlight and etc.
- G. Strenuous action of the subject or extreme electrosurgical interference may also affect the accuracy.
- H. Testee can not use enamel or other makeup.
- I. Please clean and disinfect the device after operating according to the User Manual (7.1).

6.3. Clinical restrictions

- A. As the measure is taken on the basis of arteriole pulse, substantial pulsating blood flow of subject is required. For a subject with weak pulse due to shock, low ambient/body temperature, major bleeding, or use of vascular contracting drug, the SpO₂ waveform (PLETH) will decrease. In this case, the measurement will be more sensitive to interference.
- B. For those with a substantial amount of staining dilution drug (such as methylene blue, indigo green and acid indigo blue), or carbon monoxide

hemoglobin (COHb), or methionine (Me+Hb) or thiosalicylic hemoglobin, and some with icterus problem, the SpO₂ determination by this monitor may be inaccurate.

C. The drugs like dopamine, procaine, prilocaine, lidocaine and butacaine may also be a major factor blamed for serious error of SpO₂ measure.

D. As the SpO₂ value serves as a reference value for judgement of anemic anoxia and toxic anoxia, some patients with serious anemia may also report good SpO₂ measurement.

7. Maintain, Transportation and Storage

7.1 Cleaning and Disinfecting

When using alcohol wipes to disinfect the device, please air dry or clean it with clean soft cloth.

7.2 Maintain

- Please clean and disinfect the device before using according to the User Manual (7.1).
- Please recharge the battery when the screen shows .
- Recharge the battery soon after the over-discharge. The device should be recharged every six months when it is no regular used. It can extend the battery life following this guidance.

7.3 Transportation and Storage

- The packed device can be transported by ordinary conveyance or according to transport contract. The device can not be transported mixed with toxic, harmful, corrosive material.
- The best storage environment of the device is -40 °C ~ 60 °C ambient temperature and not higher than 95% relative humidity.

8. Troubleshooting

Trouble	Possible Reason	Solution
The SpO ₂ and Pulse Rate can not be displayed normally	<ol style="list-style-type: none"> The finger is not properly positioned The patient's SpO₂ is too low to be detected 	<ol style="list-style-type: none"> Place the finger properly and try again. Try again; Go to a hospital for a diagnosis if you are sure the device works properly.
The SpO ₂ and Pulse Rate are not displayed stably	<ol style="list-style-type: none"> The finger is not placed inside deep enough. The finger is shaking or the patient is moving 	<ol style="list-style-type: none"> Place the finger properly and try again. Calm the patient down
The device can not be turned on	<ol style="list-style-type: none"> The batteries are drained or almost drained Malfunction of the device 	<ol style="list-style-type: none"> Please recharge the battery Please contact the local service center
The display suddenly turns off	<ol style="list-style-type: none"> The product will enter standby mode without 	<ol style="list-style-type: none"> Turn on the device and insert finger

	signal for 5 seconds 2. The battery is drained away or almost drained	2. Please recharge the battery
The battery can not be full charged even after 10 hours charging time.	The battery is broken	Please contact the local service center

9. Key of Symbols

Symbol	Description	Symbol	Description
	Refer to instruction manual/booklet	---	1. The finger clip falls off (no finger inserted) 2. Probe error 3. Signal inadequacy indicator
%SpO2	The pulse oxygen saturation (%)	IP22	International Protection
bpm	Pulse rate (bpm)		WEEE (2002/96/EC)
	Full-voltage		UKCA marking
	Two grid of power		One grid of power
	Low-voltage		Manufacturer
	Close the alarm sound indication		Date of manufacture
	Pause the alarm sound indication		Storage and Transport Temperature limitation
	Open the alarm sound indication		Storage and Transport Humidity limitation
	Close the pulse sound indication		Storage and Transport Atmospheric pressure limitation
	Open the pulse sound indication		This side UP
	Menu button/Exit standby mode		Fragile, handle with care
	Type BF		Keep dry
	USB		Recyclable
	Serial number		Bluetooth: ON (Bluetooth wired equipment)

R●	Record state	EC REP	European Authorised Representative
CH REP	Swiss Representative	UK REP	Authorised representative in UK
CE ₀₁₂₃	This item is compliant with Medical Device Directive 93/42/EEC of June 14, 1993, a directive of the European Economic Community.		

10. Function Specification

Information	Display Mode
The Pulse Oxygen Saturation (SpO2)	2-digit digital LCD display
Pulse Rate (bpm)	3-digit digital LCD display
Pulse Intensity (bar-graph)	Bar-graph LCD display
SpO2 Parameter Specification	
Measuring range	0%~100%, (the resolution is 1%).
Accuracy	70%~100%: ±2%, below 70% unspecified.
Average value	Calculate the Average value in every 4 measure value. The deviation between average value and true value does not exceed 1%.
Pulse Parameter Specification	
Measuring range	30bpm~250bpm, (the resolution is 1 bpm)
Accuracy	±2bpm or ±2% (select larger)
Average pulse rate	Moving calculate the Average pulse rate every 4 cardio-beat's cycle. The deviation between average value and true value does not exceed 1%
Safety Type	Interior Battery, BF Type
Pulse Intensity	
Range	Continuous bar-graph display, the higher display indicates the stronger pulse.
Battery Requirement	
Voltage 3.7 rechargeable lithium battery × 1	
Battery working life	
Charge and discharge no less than 500 times.	
Power Adapter	
Input Voltage	100 to 240 VAC, 50/60 Hz
Output voltage	DC 5 V
Output current	1000mA
Oximeter Probe	
Wavelength:660nm~905nm	
Dimensions and Weight	

Dimensions	57(L) × 32(W) × 30 (H) mm
Weight	About 50g (with the lithium battery*1)

11. Warranty and service

We provide 1 year warranty for material and manufacturing defects of the product.

The warranty does not apply:

- in case of damage caused by improper operation
- for wearing parts
- for defects that were already known to the customer at the time of purchase
- in case of the customer's own fault

The statutory warranties of the customer remain unaffected by the guarantee. For the assertion of a warranty claim within the warranty period, the customer must provide proof of purchase.

The warranty must be claimed within a period of 1 year from the date of purchase from Novidion GmbH, Fuggerstr. 30, 51149 Cologne, Germany. In the event of a guarantee claim, the customer has the right to repair the goods at our or at workshops authorised by us. Further rights are not granted to the customer on the basis of the guarantee.

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Appendix 2

Guidance and manufacture's declaration

Guidance and manufacture's declaration – electromagnetic emissions-
for all EQUIPMENT and SYSTEMS

Guidance and manufacture's declaration – electromagnetic emission		
The <i>CMS50E Pulse Oximeter</i> is intended for use in the electromagnetic environment specified below. The customer of the user of the <i>CMS50E Pulse Oximeter</i> should assure that it is used in such and environment.		
Emission test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The <i>CMS50E Pulse Oximeter</i> uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emission CISPR 11	Class B	The <i>CMS50E Pulse Oximeter</i> is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.

Guidance and manufacture's declaration – electromagnetic immunity –
for all EQUIPMENT and SYSTEMS

Guidance and manufacture's declaration – electromagnetic immunity
The <i>CMS50E Pulse Oximeter</i> is intended for use in the electromagnetic

environment specified below. The customer or the user of <i>CMS50E Pulse Oximeter</i> should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	8 kV contact 15 kV air	8 kV contact 15 kV air	Floors should be wood, concrete or ceramic tile. If floor are covered with synthetic material, the relative humidity should be at least 30%. the manufacturer may recommend the ESD precautionary procedures to user.
Power frequency (50Hz) magnetic field IEC 61000-4-8	30A/m	30A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Guidance and manufacture's declaration – electromagnetic immunity –
for EQUIPMENT and SYSTEMS that are not LIFE-SUPPORTING

Guidance and manufacture's declaration – electromagnetic immunity			
The <i>CMS50E Pulse Oximeter</i> is intended for use in the electromagnetic environment specified below. The customer or the user of <i>CMS50E Pulse Oximeter</i> should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3V(0.15MHz–80MHz),6V(in SM bands between 0.15MHz and 80MHz) 10 V/m 80 MHz to	3V(0.15MHz–80MHz),6V(in ISM bands between 0.15MHz and 80MHz) 10 V/m	Portable and mobile RF communications equipment should be used no closer to any part of the CMS50E, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = \left[\frac{3.5}{V_1} \right] \sqrt{P}$

Radiated RF EC 61000-4 3	2.7GH		$d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$ 80 MHz to 800 MHz $d = \left[\frac{7}{E_1} \right] \sqrt{P}$ 800 MHz to 2.7 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, a should be less than the compliance level in each frequency range. b Interference may occur in the vicinity of equipment marked with the following symbol: 
--------------------------------	-------	--	---

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.
 NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the CMS50E is used exceeds the applicable RF compliance level above, the CMS50E should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the CMS50E.
 b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distances between portable and mobile RF communications equipment and the EQUIPMENT or SYSTEM –

for EQUIPMENT or SYSTEM that are not LIFE-SUPPORTING

Recommended separation distances between portable and mobile RF communications equipment and the *CMS50E Pulse Oximeter*.

The *CMS50E Pulse Oximeter* is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the *CMS50E Pulse Oximeter* can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the *CMS50E Pulse Oximeter* as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter(m)		
	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.7 GHz
	$d = \left[\frac{3.5}{V_1} \right] \sqrt{P}$	$d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$	$d = \left[\frac{7}{E_1} \right] \sqrt{P}$
0.01	0.058	0.035	0.07
0.1	0.18	0.11	0.22
1	0.58	0.35	0.7
10	1.83	1.10	2.21
100	5.8	3.5	7

For transmitters rated at a maximum output power not listed above, the recommended separation distance *d* in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where *P* is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.