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## Instructions to User

Dear users, thank you very much for purchasing the Pulse Oximeter.

This Manual is written and compiled in accordance with the council directive MDD93/42/EEC for medical devices and harmonized standards. In case of modifications and software upgrades, the information contained in this document is subject to change without notice.

The Manual describes, in accordance with the Pulse Oximeter's features and requirements, main structure, functions, specifications, correct methods for transportation, installation, usage, operation, repair, maintenance and storage, etc. as well as the safety procedures to protect both the user and equipment. Refer to the respective chapters for details.

Please read the User Manual carefully before using this product. The User Manual which describes the operating procedures should be followed strictly. Failure to follow the User Manual may cause measuring abnormality, equipment damage and human injury. The manufacturer is NOT responsible for the safety, reliability and performance issues and any monitoring abnormality, human injury and equipment damage due to users' negligence of the operation instructions. The manufacturer's warranty service does not cover such faults.

Owing to the forthcoming renovation, the specific products you received may not be totally in accordance with the description of this User Manual. We apologize for this.

This product is a medical device, which can be used repeatedly.

### Warning

- ! Uncomfortable or painful feeling may appear if using the device ceaselessly, especially for the microcirculation barrier patients. It is recommended that the sensor should not be applied to the same finger for over 2 hours.
- ! For special patients, there should be a more prudent inspection in the placing process. The device can not be clipped on the edema and tender tissue.
- ! The light (infrared is invisible) emitted from the device is harmful to the eyes, so the user and maintenance personnel should not stare at the light.
- ! Testee can not use enamel or other makeup.
- ! Testee's fingernail can not be too long.
- ! Please refer to the correlative literature about the clinical restrictions and caution.
- ! This device is not intended for diagnosis.

The User Manual is published by our company. All rights reserved.

## 1. Safety

### 1.1 Instructions for Safe Operations

- ◇ Check the main unit and all accessories periodically to make sure that there is no visible damage that may affect patient's safety and monitoring performance about cables and transducers. It is recommended that the device should be inspected once a week at least. In case of obvious damage, stop using the device.
- ◇ Necessary maintenance must be performed by qualified service engineers ONLY. Users are not permitted to maintain it by themselves.
- ◇ The oximeter cannot be used together with devices not specified in User's Manual. Only accessory appointed or recommended by manufacture can be used with this device.
- ◇ This product is calibrated before leaving factory.

### 1.2 Warning

- ! Explosive hazard—DO NOT use the oximeter in environment with inflammable gas such as some ignitable anesthetic agents.
- ! DO NOT use the oximeter while the testee is being measured by MRI and CT.
- ! The person who is allergic to rubber can not use this device.
- ! The disposal of scrap instrument and its accessories and packings (including battery, plastic bags, foams and paper boxes) should follow the local laws and regulations.
- ! Please check the packing before use to make sure the device and accessories are totally in accordance with the packing list, or else the device may have the possibility of working abnormally.
- ! Please don't measure this device with functional tester for the device's related information.

### 1.3 Attention

- Keep the oximeter away from dust, vibration, corrosive substances, explosive materials, high temperature and moisture.
- If the oximeter gets wet, please stop operating it.
- When it is carried from cold environment to warm or humid environment, please do not use it immediately.
- DO NOT operate keys on front panel with sharp materials.
- High temperature or high pressure steam disinfection of the oximeter is not permitted. Refer to User Manual in the relative chapter for instructions of cleaning and disinfection.
- Do not have the oximeter immersed in liquid. When it needs cleaning, please wipe its surface with medical alcohol on soft tissue. Do not spray any liquid on the device directly.
- When cleaning the device with water, the temperature should be lower than 60°C.
- As to the fingers which are too thin or too cold, it would probably affect the normal measure of the patients' SpO<sub>2</sub> and pulse rate, please clip a thicker

finger such as thumb and middle finger deeply enough into the probe.

- Do not use the device on infant or neonatal patients.
- The product is suitable for children above four years old and adults (weight should be between 15kg and 110kg).
- The device may not work for all patients. If you are unable to achieve stable readings, discontinue use.
- The update period of data is less than 5 seconds, which is changeable according to different individual pulse rate.
- The waveform is normalized. Please read the measured value when the waveform on screen is equably and steady-going. Here this measured value is optimal value, and the waveform at the moment is the standard one.
- If some abnormal conditions appear on the screen during test process, pull out the finger and reinsert to restore normal use.
- The device has an estimated lifespan of three years from the first electrified use.
- The hanging rope attached to the device is made from Non- allergy material, if particular group are sensitive to the hanging rope, stop using it. In addition, pay attention to the use of the hanging rope, do not wear it around the neck on the purpose of avoiding harm to the patient
- This device includes a alarmfunction, users can check this according to chapter 7.
- The device has the function of limits alarming, when the measured data is beyond the highest or lowest limit, the device would start alarming automatically on the premise of the alarming function is on.
- The device has the function of alarming, this function can either be paused, or closed (default setting) for good. This function could be turned on through menu operation if you need. Please check the chapter 7 as a reference.
- The instrument does not have low-voltage alarm function, it only shows the low-voltage. Please change the battery when the battery energy is used out.
- Batteries must be removed if the device is going to be stored for more than one month, or else batteries may leak.
- A flexible circuit connects the two parts of the device. Do not twist or pull on the connection.

#### **1.4 Indication for use**

The Fingertip Pulse Oximeter is a non-invasive device intended for the spot-check of oxygen saturation of arterial hemoglobin (SpO<sub>2</sub>) and the pulse rate of adult and pediatric patients in home and hospital environments (including clinical use in internist/surgery, anesthesia, intensive care etc.). This device is not intended for continuous monitoring.

## **2. Overview**

The pulse oxygen saturation is the percentage of HbO<sub>2</sub> in the total Hb in the blood, in other words the O<sub>2</sub> concentration in blood. It is an important bio-parameter for respiration. A number of diseases relating to the respiratory system

may cause a decrease of SpO<sub>2</sub> in the blood. Serious symptoms might bring danger to person's life. Therefore, prompt information of patients' SpO<sub>2</sub> is of great help for the doctor to discover the potential danger, and is of great importance in the clinical medical field.

The Pulse Oximeter features in small volume, low power consumption, convenient operation and being portable. It is only necessary for patients to put one of his fingers into a probe for diagnosis, and a display screen will directly show the measured value of pulse oxygen saturation.

### 2.1 Classification

Class IIb, (MDD93/42/EEC IX Rule 10)

### 2.2 Features

- A. Operation of the product is simple and convenient.
- B. The product is small in volume, light in weight and convenient for carrying.
- C. Low power consumption
- D. The product will automatically be powered off when there is no signal for 5 seconds.

### 2.3 Major Applications and Scope of Application

The Pulse Oximeter can be used in measuring the pulse oxygen saturation and pulse rate through finger. The product is suitable for being used in family, hospital, oxygen bar, community healthcare, physical care in sports (It can be used before or after doing sports and it is not recommended to use the device during the process of exercises) and more.

The product is not suitable for use in continuous supervision for patients.

The problem of overrating would emerge when the patient is suffering from toxicosis caused by carbon monoxide, the device is not recommended to be used under this circumstance.

### 2.4 Environment Requirements

Storage Environment

- a. Temperature : -40°C~+60°C
- b. Relative humidity : ≤95%
- c. Atmospheric pressure : 500hPa~1060hPa

Operating Environment

- a. Temperature: 10°C ~40°C
- b. Relative Humidity : ≤75%
- c. Atmospheric pressure: 700hPa~1060hPa

## 3. Principle

### 3.1 Principle of Measurement

The principle of the Oximeter is as follows: An experience formula of data process is established taking use of Lambert Beer Law according to Spectrum Absorption Characteristics of Reductive Hemoglobin (Hb) and Oxyhemoglobin (HbO<sub>2</sub>) in glow & near-infrared zones. The operation principle of the device is: Photoelectric Oxyhemoglobin Inspection Technology is adopted in accordance with Capacity Pulse Scanning & Recording Technology, so that two beams of different

wavelengths of lights can be focused onto a human nail tip through perspective clamp finger-type sensor. The measured signal can be obtained by a photosensitive element, information acquired through which will be shown on screen through treatment in electronic circuits and microprocessor.


Figure 1 Operating principle

### 3.2. Caution

1. The finger should be placed properly (see the attached illustration of this manual, Figure 4), or else it may cause inaccurate measurement.
2. The SpO<sub>2</sub> probe and photoelectric receiving tube should be arranged in a way with the subject's arteriole in a position there between
3. The SpO<sub>2</sub> probe should not be used at a location or limb tied with arterial canal or blood pressure cuff or receiving intravenous injection.
4. Make sure the optical path is free from any optical obstacles like rubberized fabric.
5. Excessive ambient light may affect the measuring result. It includes fluorescent lamps, dual ruby light, infrared heater, direct sunlight and etc.
6. Strenuous action of the subject or extreme electro-surgical interference may also affect the accuracy.
7. Testee can not use enamel or other makeup.

### 3.3. Clinical restrictions

**A.** As the measurements are taken on the basis of arteriole pulse, substantial pulsating blood flow of subject is required. For a subject with weak pulse due to shock, low ambient/body temperature, major bleeding, or use of vascular contracting drug, the SpO<sub>2</sub> waveform (PLETH) will decrease. In this case the measurement will be more sensitive to interference.

**B.** For those with a substantial amount of staining dilution drug (such as methylene blue, indigo green and acid indigo blue), or carbon monoxide hemoglobin (COHb), or methionine (Me+Hb) or thiosalicylic hemoglobin, and some with icterus problem, the SpO<sub>2</sub> determination by this monitor may be inaccurate.

**C.** The drugs like dopamine, procaine, prilocaine, lidocaine and butacaine may also be a major factor blamed for serious error of SpO<sub>2</sub> measure.

**D.** As the SpO<sub>2</sub> value serves as a reference value for judgement of anemic anoxia and toxic anoxia, some patients with serious anemia may also report good SpO<sub>2</sub> measurement.

## 4. Technical Specifications

- A. Display Format:** TFT Display  
**SpO<sub>2</sub> Measuring range:** 0~100 %

**Pulse** Measuring range: 30bpm ~250bpm

**Perfusion Index Measuring Range:** 0% ~ 20%

**Pulse Wave Display:** columniation display and the waveform display

- B. Power Requirements:** 2 ×1.5V AAA alkaline battery (or using the rechargeable battery instead), adaptable range: 2.6V~3.6V.
- C. Power Consumption:** Smaller than 30mA
- D. Resolution:** SpO<sub>2</sub>: 1%, Pulse: 1bpm, Perfusion Index: 0.1%.
- E. Measurement Accuracy:** ±2% in stage of 70%-100% SpO<sub>2</sub>, and meaningless when stage being smaller than 70%. ± 2 bpm during the pulse rate range of 30-99 bpm and ± 2% during the pulse rate range of 100~250 bpm
- F. Measurement Performance in weak filling condition:** SpO<sub>2</sub> and pulse rate can be shown correctly when pulse-filling ratio is 0.4%. SpO<sub>2</sub> error is ±4%, pulse rate error is ±2 bpm or ±2% (select larger).
- G. Resistance to surrounding light:** The deviation between the value measured in the condition of man-made light or indoor natural light and that of darkroom is less than ±1%.
- H.** It is equipped with a function switch. The Oximeter can be powered off in case no finger in the Oximeter within 5 seconds.
- I. Optical Sensor**  
 Red light (wavelength is 660nm, 6.65mW)  
 Infrared (wavelength is 880nm,6.75mW)
- J. Adjustable alarm range**  
 SpO<sub>2</sub> : 0%~100%  
 Pulse Rate : 0bpm~254bpm

## 5. Accessories

Hanging rope, pouch, silicone protector, batteries (optional), user manual

## 6. Installation

### 6.1 View of the front panel


Fig. 2 Front view

Fig. 3 Batteries installation

### 6.2 Battery

Step 1. Refer to Figure 3 and insert the two AAA size batteries properly in the right direction.

Step 2. Replace the cover.

Please take care when you insert the batteries for the improper insertion

may damage the device.

### 6.3 Mounting the Hanging Rope

Step 1. Put the end of the rope through the hole.

Step 2. Put another end of the rope through the first one and then tighten it.

## 7. Operating Guide

### A. Measurement

- 1) Insert the two batteries properly in the correct direction and then replace the cover.
- 2) Open the clip as shown in Figure 4.

Figure 4 Finger Position

- 3) Put the patient's finger into the rubber cushions of the clip (make sure the finger is in the right position), and then clip the finger.
- 4) Press the switch button once on front panel.
- 5) Do not shake the finger and keep the patient at ease during the process. Meanwhile, human body is not recommended in movement status.
- 6) Get the information directly from screen display.
- 7) The button  ·  ·  has three functions. When the device is off, press the button to turn it on. When the device is on, press the button shortly to enter the menu. In the state of alarm, press the button long to pause the alarm for 60s.

### B. Change display direction:

The display direction of the screen changes automatically.

### C. Pause sound prompt

1. Sound prompt includes over-limit prompt or finger out prompt.
2. In sound prompt "ON" state, when the sound prompt occurs, short press the button to pause the sound prompt, and it will resume automatically after about 60s.
3. If you want to turn off the sound prompt permanently, please set it in menu.

### D. Menu operation

#### a. Main menu

 Pulse	off	→ Pulse Sound on / off
Alarm	off	→ Alarm on / off
Alarm Set		→ Sound prompt menu
Exit		→ Exit main menu

**Pulse:** When you enter the main menu, „Pulse“ is already selected. Press the button for a few seconds to turn on the pulse sound. Long press again to switch off the pulse sound. Press the button long to turn the pulse sound on. Long press again to switch off the pulse sound.

**Alarm:** Press the button briefly to select „Alarm“. Press the button for a few seconds to turn on the alarm prompt on. Long press again to switch off the alarm prompt.

**Alarm Set:** Press the button briefly to select „Alarm Set“. Press the button for a few seconds to go into the sound prompt menu. For details on setting the limit values, see c. Set limit values.

**Exit:** Short press the button to select “Exit”, then long press the button to exit the main menu and return to the measurement interface.

#### b. Sound prompt menu

- Set Pulse Limits
- Set SpO2 Limits
- Exit Sound Prompt menu

**Pulse Alarm:** Press the button briefly to select „Pulse Alarm“. Press the button for a few seconds to go to the Pulse limit value settings, see c. Set limit values for details.

**SpO<sub>2</sub> Alarm:** Press the button briefly to select „ SpO<sub>2</sub> Alarm“. Press the button for a few seconds to go to the SpO2 limit value settings, see c. Set limit values for details.

**Exit:** Briefly press the button to select “Exit”, then long press the button to exit the alarm menu and return to the main menu.

#### c. Set limit values

- Direction: up, down
- Upper limit prompt
- Lower limit prompt

**Dir:** When you enter the limit value settings, „Dir“ is already selected. Long press the button to switch between up (increase the value) and down (decrease the value).

**High:** Briefly press the button to select “High”. Long press the button to set the upper limit prompt of pulse or SpO<sub>2</sub>. Make sure that the direction in which you want to change the value is correctly set to up or down. Lower limit can not exceed the upper limit.

**Low:** Briefly press the button to select “Low”. Long press the button to set the lower limit prompt of pulse or SpO<sub>2</sub>. Make sure that the direction in which you want to change the value is correctly set to up or down. Lower limit can not exceed the upper limit.

**Exit:** Briefly press the button to select “Exit”, then long press the button to exit and return to the alarm menu.

## 8. Cleaning, Transportation and Storage

- Please change the batteries when the low-voltage is being displayed on the screen.
- Please clean the surface of the device before using. Wipe the device with

medical alcohol first, and then let it dry on air or clean it with dry clean fabric.




- Using the medical alcohol to disinfect the product after use, prevent from cross infection for next time use.
- Please take out the batteries if the oximeter is not in use for a long time.
- The best storage environment of the device is  $-40^{\circ}\text{C}$  to  $60^{\circ}\text{C}$  ambient temperature and not higher than 95% relative humidity.

-  High-pressure sterilization cannot be used on the device.
-  Do not immerse the device in liquid.
-  It is recommended that the device should be kept in a dry environment. Humidity may reduce the useful life of the device, or even damage it.

## 9. Troubleshooting

Trouble	Possible Reason	Solution
The SpO2 and Pulse Rate can not be displayed normally	<ol style="list-style-type: none"> <li>1. The finger is not properly positioned</li> <li>2. The patient's SpO2 is too low to be detected.</li> </ol>	<ol style="list-style-type: none"> <li>1. Place the finger properly and try again.</li> <li>2. Try again; Go to a hospital for a diagnosis if you are sure the device works properly.</li> </ol>
The SpO2 and Pulse Rate are not displayed stably	<ol style="list-style-type: none"> <li>1. The finger is not placed inside deep enough.</li> <li>2. The finger is shaking or the patient is moving.</li> </ol>	<ol style="list-style-type: none"> <li>1. Place the finger properly and try again.</li> <li>2. Calm thw patient down</li> </ol>
The device can not be turned on	<ol style="list-style-type: none"> <li>1. The batteries are drained or almost drained.</li> <li>2. The batteries are not inserted properly.</li> <li>3. Malfunction of the device.</li> </ol>	<ol style="list-style-type: none"> <li>1. Change batteries</li> <li>2. Reinstall batteries.</li> <li>3. Please contact the local service center.</li> </ol>
The display suddenly turns off	<ol style="list-style-type: none"> <li>1. The device will power off automatically without signal within 5 seconds.</li> <li>2. The battery is almost drained.</li> </ol>	<ol style="list-style-type: none"> <li>1. Turn on the device and insert finger.</li> <li>2. Change batteries</li> </ol>

## 10. Key of Symbols

Symbol	Description	Symbol	Description
	Type BF applied part		Menu/ Power button
	Refer to instruction manual/booklet		Sound prompt on
%SpO2	Pulse oxygen saturation (%)		Pulse sound on
bpm	Pulse rate (bpm)		Serial number
PI	Perfusion Index (%)		WEEE (2002/96/EC)
	Low battery		UKCA marking
	<ol style="list-style-type: none"> <li>1. No finger inserted</li> <li>2. An indicator of signal inadequacy</li> </ol>	IP22	International Protection
	Battery anode / Battery cathode		Manufacturer
	Date of manufacture		European authorised representative
	UK Authorised Representative		Swiss authorised

	This item is compliant with Medical Device Directive 93/42/EEC of June 14, 1993, a directive of the European Economic Community.	representative
--	--	----------------

## 11. Function Specification

Information		Display Mode
The Pulse Oxygen Saturation ( SpO <sub>2</sub> )		LCD
Pulse Rate ( PR )		LCD
Pulse Intensity (bar-graph)		LCD
Pulse wave		LCD
SpO <sub>2</sub> Parameter Specification		
Measuring range	0%~100%, (the resolution is 1%).	
Accuracy	70%~100%:±2% ,Below 70% unspecified.	
Optical Sensor	Red light (wavelength is 660nm) Infrared (wavelength is 880nm)	
Pulse Parameter Specification		
Measuring range	30bpm~250bpm, (the resolution is 1bpm)	
Accuracy	±2bpm or±2% (select larger)	
Pulse Intensity		
Range	Continuous bar-graph display, the higher display indicates the stronger pulse.	
Battery Requirement		
1.5V (AAA size) alkaline batteries × 2		
Battery useful life		
Two batteries can work continuously for 20 hours		
Dimensions and Weight		
Dimensions	57(L) × 31(W) × 32 (H) mm	
Weight	About 50g (with batteries)	

## 12. Warranty and service

We provide 1 year warranty for material and manufacturing defects of the product.

The warranty does not apply:

- in case of damage caused by improper operation
- for wearing parts
- for defects that were already known to the customer at the time of purchase
- in case of the customer's own fault

The statutory warranties of the customer remain unaffected by the guarantee.

For the assertion of a warranty claim within the warranty period, the customer must provide proof of purchase.

The warranty must be claimed within one year from the date of purchase to:

Novidion GmbH, Fuggerstr. 30, 51149 Cologne, Germany.

In case of a warranty claim, the customer has the right to have the goods repaired by us or in workshops authorized by us.

The customer is not granted any further rights based on the warranty.

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### **CONTEC MEDICAL SYSTEMS CO, LTD**

ADD: No.112 Qinhuang West Street, Economic & Technical Development Zone,  
Qinhuangdao, Hebei Province, PEOPLE'S REPUBLIC OF CHINA

### **Prolinx GmbH**

ADD: Brehmstr. 56, 40239 Duesseldorf, Germany  
Tel: 0049 211 3105 4698  
E-mail: med@eulinx.eu

### **MedPath Limited**

ADD: 27 Old Gloucester Street, London,  
United Kingdom, WC1N 3AX  
E-mail: UK@medpath.pro

### **Tinovamed GmbH**

ADD: Zielmattenring 6, 4563 Gerlafingen, Switzerland  
E-mail: info@tinovamed.ch

**Bei Fragen wenden Sie sich bitte an / If you have any questions, please contact /**  
Si tiene alguna pregunta, póngase en contacto con nosotros / Si vous avez une  
question s'il vous plaît contacter/Se si dispone di una domanda si prega di  
contattare / Als u vragen heeft, kunt u contact opnemen met / W razie  
jakichkolwiek pytań prosimy o kontakt / Om du har några frågor, vänligen  
kontakta:

Novidion GmbH, Fuggerstr. 30, 51149 Köln, Deutschland  
Tel.: +49 2203 – 9885 200, Fax: +49 2203 – 9885 206  
info@novidion.de, www.pulox.de  
WEEE: DE24355330

**Sollten Sie beim Lesen dieser Anleitung Probleme haben, so können Sie diese  
auch auf unserer Website downloaden / In case you have any problems reading  
this manual, you can download it from our website / En caso de que tenga algún  
problema para leer este manual, puede descargarlo de nuestro sitio web / Si vous  
rencontrez des problèmes pour lire ce manuel, vous pouvez le télécharger depuis  
notre site Web / In caso di problemi nella lettura di questo manuale, è possibile  
scaricarlo dal nostro sito Web / Als u problemen hebt met het lezen van deze  
instructies, kunt u ze ook downloaden van onze website / Jeśli masz problemy z  
odczytaniem tej instrukcji, możesz ją również pobrać z naszej strony internetowej  
/ Om du har problem med att läsa dessa instruktioner kan du också ladda ner dem  
från vår webbplats:**

## Appendix

### Guidance and manufacture's declaration-electromagnetic emission for all EQUIPMENT and SYSTEMS

Guidance and manufacture's declaration –electromagnetic emission			
The <i>Pulse Oximeter</i> is tended for use in the electromagnetic environment specified below. The customer of the user of the <i>Pulse Oximeter</i> should assure that it is used in such an environment.			
Emission test	compliance	Electromagnetic environment-guidance	
RF emissions CISPR 11	Group 1	The <i>Pulse Oximeter</i> uses RF energy only for their internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.	
RF emissions CISPR 11	Class B	The <i>Pulse Oximeter</i> is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.	
Harmonic emissions IEC 61000-3-2	Not applicable		
Voltage fluctuations/ flicker emission IEC 61000-3-3	Not applicable		

### Guidance and manufacture's declaration-electromagnetic immunity for all EQUIPMENT and SYSTEMS

Guidance and manufacture's declaration-electromagnetic immunity			
The <i>Pulse Oximeter</i> is intended for use in the electromagnetic environment specified below. The the user of <i>Pulse Oximeter</i> should assure that it is used in such an environment.			
Immunity test	IEC60601 test level	Compliance level	Electromagnetic environment-guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6KV contact ±8KV air	±6KV contact ±8KV air	Floors should be wood, concrete or ceramic tile. If floor are covered with synthetic material, the relative humidity should be at least 30%.
Power frequency (50Hz) magnetic field IEC 61000-4-8	3A/m	3A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment

### Guidance and manufacture's declaration-electromagnetic immunity for EQUIPMENT and SYSTEMS that are not LIFE-SUPPORTING

**Guidance and manufacture's declaration-electromagnetic immunity**

The *Pulse Oximeter* is intended for use in the electromagnetic environment specified below. The customer or the user of *Pulse Oximeter* should assure that it is used in such an environment.

Immunity test	IEC60601 test level	Compliance level	Electromagnetic environment - guidance
Radiated RF ICE 61000-4-3	3V/m 80MHz to 2.5GHz	3V/m	<p>Portable and mobile RF communication equipment should be used no closer to any part of the <i>Pulse Oximeter</i>, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>recommended separation distance</p> <p>80MHz to 800MHz</p> $d = \left[ \frac{3.5}{E_1} \right] \sqrt{P}$ <p>800MHz to 2.5GHz</p> $d = \left[ \frac{7}{E_1} \right] \sqrt{P}$ <p>Where <math>P</math> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and <math>d</math> is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,<sup>a</sup> should be less than the compliance level in each frequency range.<sup>b</sup> Interference may occur in the vicinity of equipment marked with the following symbol: </p>

NOTE 1 At 80MHz and 800MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

<sup>a</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the

measured field strength in the location in which The Pulse Oximeter is used exceeds the applicable RF compliance level above, the Pulse Oximeter should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Pulse Oximeter.

<sup>b</sup> Over the frequency range 150 KHz to 80 MHz, field strengths should be less than 3V/m.

**Recommended separation distances between portable and mobile RF communications equipment and the EQUIPMENT or SYSTEM for EQUIPMENT or SYSTEM that not LIFE-SUPPORTING**

**Recommended separation distances between portable and mobile RF communications equipment and the Pulse Oximeter**

The *Pulse Oximeter* is intended for use in the electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the *Pulse Oximeter* can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the *Pulse Oximeter* as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150KHz to 80MHz $d = \left[ \frac{3.5}{E_1} \right] \sqrt{P}$	80MHz to 800MHz $d = \left[ \frac{3.5}{E_1} \right] \sqrt{P}$	800MHz to 2.5GHz $d = \left[ \frac{7}{E_1} \right] \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.74
1	1.17	1.17	2.33
10	3.69	3.69	7.38
100	11.67	11.67	23.33

For transmitters rated at a maximum output power not listed above, the recommended separation distance in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80MHz and 800MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.