# 3M Tegaderm™ CHG

(en)	Chlorhexidine Gluconate Gel Pad .....	1
(fr)	Tampon de gel au gluconate de chlorhexidine .....	4
(de)	Chlorhexidingluconat-Gelkissen .....	7
(it)	Pad in gel con clorexidina gluconato .....	10
(es)	Parche de gel con gluconato de clorhexidina.....	13
(nl)	Gelpad, chloorhexidinegluconaat.....	16
(sv)	Klohexidinglukonat Gel Pad.....	19
(da)	Gel Pad med klorhexidin-glukonat .....	22
(no)	Gelpute med klorheksidinglukonat.....	25
(fi)	Klooriheksidiiniglukonaattia sisältävä geelityyny.....	28
(pt)	Almofada de Gel de Gluconato de Cloroheixidina .....	31
(el)	Επίθεμα Γέλης Γλυκονικής Χλωρεξιδίνης .....	34
(pl)	Podkładka żelowa z glukonianem chlorheksydyny .....	37
(hu)	Klórhexidin-glükonátos gélpárna.....	40
(cs)	Gelová podložka z chlorhexidin-glukonátu .....	43
(sk)	Gélový vankúšik z chlórhexidínglukonátu .....	46
(sl)	Blazinica z gelom in klorheksidin glukonatom .....	49
(et)	Kloorheksidiinglükonaadiga geelpadjand .....	52
(lv)	Hlorheksidīna glikonāta gela pārklājs.....	55
(lt)	Chlorheksidino gliukonato gelio įdėklas .....	58
(ro)	Tampon cu gel cu gluconat de clorhexidină.....	61
(uk)	Гелева подушечка з хлоргексидину глюконатом .....	64
(hr)	Gelasti jastučić od klorheksidin glukonata .....	67
(bg)	Гел подложка с хлорхексидин глюконат .....	70
(sr)	Gelasti jastučić od klorheksidin glukonata .....	73
(tr)	Klorheksidin Glukonat Jel Ped .....	76
(kk)	Хлоргексидин глюконаты негізіндегі гель жастықша .....	79
(th)	แผ่นเจลคลอรัลเฮกซิดีนกลูโคเนท .....	83
(sq)	Fletë me xhel me glukonat klorheksidine.....	86
(mk)	Гел влошка со хлорхексидин глюконат .....	89

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3M and Tegaderm are trademarks of 3M.  
Used under license in Canada.  
© 2023, 3M. Unauthorized use prohibited.  
All rights reserved.

Issue Date: 2023-08  
34-8730-0048-2

3M Deutschland GmbH  
Health Care Business  
Carl-Schurz-Str. 1  
41453 Neuss, Germany

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**Chlorhexidine Gluconate Gel Pad****Product Description**

3M™ Tegaderm™ CHG Chlorhexidine Gluconate Gel Pad is used to protect vascular access sites. The CHG Gel Pad is breathable and transparent, allowing continuous site observation.

The antimicrobial CHG Gel Pad consists of an integrated gel pad containing 2% w/w chlorhexidine gluconate (CHG), a well-known antiseptic agent with broad spectrum antimicrobial and antifungal activity. The CHG Gel Pad absorbs fluid.

*In vitro* testing (time kill and zone of inhibition) demonstrates that the CHG Gel Pad has an antimicrobial effect against a variety of gram-positive and gram-negative bacteria and yeast.

**Intended Use/Intended Purpose**

3M™ Tegaderm™ CHG Chlorhexidine Gluconate Gel Pad can be used to protect catheter sites. Common applications include protecting intravascular catheters and percutaneous devices. Tegaderm™ CHG is intended to reduce skin colonization and catheter colonization and to suppress regrowth of microorganisms commonly related to bloodstream infections. Tegaderm™ CHG is intended to reduce catheter-related bloodstream infections (CRBSI) in patients with central venous or arterial catheters.

**Warnings**

DO NOT PUNCTURE CHG GEL WITH NONCORING NEEDLE.

DO NOT USE TEGADERM™ CHG DRESSINGS ON PREMATURE INFANTS OR INFANTS YOUNGER THAN 2 MONTHS OF AGE. USE OF THIS PRODUCT ON PREMATURE INFANTS MAY RESULT IN HYPERSENSITIVITY REACTIONS OR NECROSIS OF THE SKIN.

THE SAFETY AND EFFECTIVENESS OF TEGADERM™ CHG DRESSINGS HAS NOT BEEN EVALUATED IN CHILDREN UNDER 18 YEARS OF AGE.

FOR EXTERNAL USE ONLY. DO NOT ALLOW THIS PRODUCT TO CONTACT EARS, EYES, MOUTH OR MUCOUS MEMBRANES.

DO NOT USE THIS PRODUCT ON PATIENTS WITH KNOWN HYPERSENSITIVITY TO CHLORHEXIDINE GLUCONATE. THE USE OF CHLORHEXIDINE GLUCONATE CONTAINING PRODUCTS HAS BEEN REPORTED TO CAUSE IRRITATIONS, SENSITIZATION, AND GENERALIZED ALLERGIC REACTIONS. IF ALLERGIC REACTIONS OCCUR, DISCONTINUE USE IMMEDIATELY, AND IF SEVERE, CONTACT A PHYSICIAN.

Hypersensitivity reactions associated with topical use of chlorhexidine gluconate have been reported in several countries. The most serious reactions (including anaphylaxis) have occurred in patients treated with lubricants containing chlorhexidine gluconate, which were used during urinary tract procedures. Caution should be used when using chlorhexidine gluconate containing preparations, and the patient should be observed for the possibility of hypersensitivity reactions.

**Precautions**

3M™ Tegaderm™ CHG Chlorhexidine Gluconate Gel Pad should not be placed over infected wounds. It is not intended to be used as a treatment of percutaneous device-related infections.

In the case of clinical wound infection, systemic antibacterials should be used if indicated.

Any active bleeding at the insertion site should be stabilized before applying the CHG Gel Pad.

Do not stretch the CHG Gel Pad during application. Mechanical skin trauma may result if the CHG Gel Pad is applied with tension. The skin should be clean, dry, and free of detergent residue to prevent skin irritation and to ensure good adhesion.

Do not reuse. Reuse may result in compromising product integrity and lead to product failure.

Use care and cover the CHG Gel Pad to protect from water when showering.

The CHG Gel Pad has not been tested for use during radiation therapy.

**Clinical Trial Results:** A randomized, controlled clinical trial consisting of 1879 subjects with 4163 central venous and arterial catheter insertion sites was conducted at 11 hospitals (1). Results showed that the use of Tegaderm™ CHG resulted in a statistically significant 60% reduction in the incidence of catheter-related bloodstream infections ( $P=0.02$ ). Study results also demonstrate a statistically significant reduction in skin colonization ( $P<0.001$ ) and catheter colonization ( $P<0.0001$ ) in the chlorhexidine vs. non-chlorhexidine group.

Variable	Non-chlorhexidine vs. chlorhexidine dressings (941 patients/2055 catheters vs. 938 patients/2108 catheters)
<b>Catheter-related bloodstream infection</b>	
Incidence densities (n per 1000 catheter-days)	1.3 vs. 0.5
Hazard ratio	0.402 [0.186 to 0.868], $P=0.02$
<b>Catheter colonization</b>	
Incidence densities (n per 1000 catheter-days)	10.9 vs. 4.3
Hazard ratio	0.412 [0.306 to 0.556], $P<0.0001$

(1) Timsit JF et al (Randomized Controlled Trial of Chlorhexidine Dressing and Highly Adhesive Dressing for Preventing Catheter-Related Infections in Critically Ill Adults  
*Am. J. Respir. Crit. Care Med.* 2012; 186:1272-1278

A summary of safety and clinical performance (SSCP) is available in the European database on medical devices (EUDAMED), where it is linked to the basic UDI-DI. The EUDAMED public website is <https://ec.europa.eu/tools/eudamed>. The CHG Gel Pad's basic UDI-DI is 0608223840101000000039AM.


**Instructions for Use:** Refer to packaging figures and those included in package insert. Failure to follow the manufacturer's instructions for use may result in complications including skin irritation and/or maceration. If the sterile packaging is damaged or unintentionally opened before use, discard the product, do not use. This product is to be used by healthcare providers.


**Site Preparation:** Prepare site according to institution protocol. Clipping of hair at the site may improve CHG Gel Pad adhesion. Shaving is not recommended. The skin should be clean, dry and free of detergent residue.

Any active bleeding at the insertion site should be stabilized before application of CHG Gel Pad.

**Instructions for Use**

1. Open package and remove sterile CHG Gel Pad using aseptic technique.
2. Prepare skin surrounding the percutaneous device per facility protocol allowing all preps to dry completely before application to prevent skin irritation and ensure good adhesion.
3. Remove liners according to figures then apply CHG Gel Pad to the device site using firm pressure to enhance adhesion.
4. To prepare for an easy removal when used with an additional dressing, place CHG Gel Pad around the device ensuring that the device rests upon the slit of the gel pad.
5. Per facility protocol, secure the device and CHG Gel Pad to the skin with a breathable transparent semipermeable dressing.

### Site Care

1. The site should be observed daily for signs of infection or other complications. If infection is suspected, remove the CHG Gel Pad, inspect the site directly, and determine appropriate medical intervention. Infection may be signaled by fever, pain, redness, swelling, or unusual odor or discharge.
2. Inspect the CHG Gel Pad daily and change as necessary, in accordance with facility protocol; at a minimum of every 7 days, per Centers for Disease Control and Prevention (CDC) recommendations. CHG Gel Pad changes may be needed more frequently with highly exudative sites.
3. Use care and cover the CHG Gel Pad to protect from water when showering.

The 3M™ Tegaderm™ CHG Chlorhexidine Gluconate Gel Pad should be changed as necessary:

- If the CHG Gel Pad or dressing becomes loose, soiled or compromised in any way
- If the site is obscured or no longer visible
- If there is visible drainage outside the CHG Gel Pad
- If the CHG Gel Pad appears to be saturated or overly swollen\*

\*Note: To test if the CHG Gel Pad is fully saturated, lightly press down on a corner of the gel pad with your finger. If the CHG Gel Pad remains displaced once your finger is removed, the CHG Gel Pad should be changed.

\*Note: Tegaderm™ CHG Gel Pad is not intended to absorb large quantities of blood or fluid.

### Removal

1. Remove transparent semipermeable dressing in a low and slow manner, stabilizing the intravascular catheter or percutaneous device.
2. CHG Gel Pad may remain attached to the transparent semipermeable dressing during removal. Stretch release may be required to help release film dressing from the CHG Gel Pad.
3. If needed, use sterile alcohol swabs or wipes, or sterile solutions (i.e., sterile water or normal saline) to facilitate removal of the CHG Gel Pad.
4. Dispose of CHG Gel Pad per facility protocol.

### Storage/Shelf Life

For best results, store in a cool, dry place. For shelf life, refer to the expiration date on the package. Sterility of the product is guaranteed unless individual package is damaged or open.

### How Supplied/Ordering Information

Catalog #	Size	Average Amount of CHG per CHG Gel Pad (mg based on gel pad size)
1664R 	Dimensions 6.2 cm x 4.9 cm (2 7/16 in x 1 15/16 in)	30

Please report a serious incident that has occurred in relation to the product to 3M and the local competent authority.

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## Symbol Glossary

Symbol Title	Symbol	Description and Reference
Manufacturer		Indicates the medical device manufacturer. Source: ISO 15223, 5.1.1
Date of manufacture		Indicates the date when the medical device was manufactured. Source: ISO 15223, 5.1.3
Use-by date		Indicates the date after which the medical device is not to be used. Source: ISO 15223, 5.1.4
Batch code		Indicates the manufacturer's batch code so that the batch or lot can be identified. Source: ISO 15223, 5.1.5
Catalogue number		Indicates the manufacturer's catalogue number so that the medical device can be identified. Source: ISO 15223, 5.1.6
Sterilized using ethylene oxide		Indicates a medical device that has been sterilized using ethylene oxide. Source: ISO 15223, 5.2.3
Do not resterilize		Indicates a medical device that is not to be resterilized. Source: ISO 15223, 5.2.6
Do not use if package is damaged and consult instructions for use		Indicates a medical device that should not be used if the package has been damaged or opened and that the user should consult the instructions for use for additional information. Source: ISO 15223, 5.2.8
Temperature limit		Indicates the temperature limits to which the medical device can be safely exposed. Source: ISO 15223, 5.3.7
Do not re-use		Indicates a medical device that is intended for one single use only. Source: ISO 15223, 5.4.2
Consult instructions for use or consult electronic instructions for use		Indicates the need for the user to consult the instructions for use. Source: ISO 15223, 5.4.3

Does not contain natural rubber latex		Indicates there is no presence of dry natural rubber nor natural rubber latex as a material of construction within the medical device or the packaging of a medical device. Source: ISO 15223, 5.4.5 and Annex B
Medical Device		Indicates the item is a medical device. Source: ISO 15223, 5.7.7
Contains a medicinal substance		Indicates a medical device that contains or incorporates a medicinal substance. Source: ISO 15223, 5.4.7
Unique device identifier		Indicates a carrier that contains unique device identifier information. Source: ISO 15223, 5.7.10
Importer		Indicates the entity importing the medical device into the locale. Source: ISO 15223, 5.1.8
Single sterile barrier system		Indicates a single sterile barrier system. Source: ISO 15223, 5.2.11
UKCA Mark		Indicates conformity to all applicable regulations and/or directives in the United Kingdom (UK), for products placed on the market in Great Britain (GB) with UK Approved Body involvement.
Swiss Authorised Representative		Indicates the authorised representative in Switzerland. Source: Swissmedic.ch
CE Mark		Indicates conformity to all applicable European Union Regulations and Directives with notified body involvement.
Green Dot		Indicates a financial contribution to national packaging recovery company per European Directive No. 94/62 and corresponding national law. Packaging Recovery Organization Europe.

For more information, see [HCBG regulatory.3M.com](https://www.hcbg.gov.sg/regulatory/3M.com)

(en) Manufacturer  
 (fr) Fabricant  
 (de) Hersteller  
 (it) Produttore  
 (es) Fabricante  
 (nl) Fabrikant  
 (sv) Tillverkare  
 (da) Producent  
 (no) Produsent  
 (fi) Valmistaja  
 (pt) Fabricante  
 (el) Κατασκευαστής  
 (pl) Wytwórca  
 (hu) Gyártó  
 (cs) Výrobce  
 (sk) Výrobca  
 (sl) Proizvajalec  
 (et) Tootja  
 (lv) Ražotājs  
 (lt) Gamintojas  
 (ro) Producător  
 (uk) Виробник  
 (hr) Proizvođač  
 (bg) Производител  
 (sr) Proizvođač  
 (tr) Üretici  
 (kk) Өндіруші  
 (th) ผู้ผลิต  
 (sq) Prodhuesi  
 (mk) Производител

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