- Carefully position the two lowest hook straps over the leg and fasten them without applying compression. The hook fasteners should overlap as little as possible but do not leave a gap between the hook fasteners. Press the ends of the fasteners firmly to the sides. Continue positioning the remaining hook fasteners in the same way and fix them in place.
- Apply the Accutabs and complete applying the hook fasteners according to the instructions in the package leaflet. Your prescribing physician should tell you which compression range is best for you.

Doffing the Compreflex® Standard Calf

For doffing, loosen each hook fastener, roll their ends backwards, and secure them there so that they cannot get tangled. Remove the product from the leg and take off the liner.

Storage, wearing period, and service life

Unless otherwise prescribed by your physician, the compression product can be worn during the day and overnight. However, remove the product once a day to clean the skin and care for any wounds. The product should then be reapplied according to the instructions for use. How often the product needs to be removed for skin checks, hygiene, and/or wound care may vary according to the individual needs of patients and the prescriber's recommendations. A liner should be worn under every compression product. Wearing and washing compression products every day may cause the medical pressure gradient and elasticity of your wrap to diminish. The shelf life date is indicated on the packaging with an hourglass symbol. We guarantee maximum effect for the subsequent six months of wearing.

Care and storage

Follow the instructions for care on the packaging and/or the textile label. Close all hook fasteners before washing. Wash the product in a mesh laundry bag on gentle cycle using a mild detergent. Do not use bleach or fabric softener. Dry the product on low/no heat in a dryer or

air dry (air drying is recommended to extend wrap life). Store the product at room temperature and in a dry place away from light. Keep away from direct heat.
Wash liners before first use.

Disposal

There are no special requirements for disposing of the medical compression product. They can be disposed of as residual or household waste. Please comply with local regulations for the disposal of textiles.

Material abbreviations

37% Polyamide
27% Polyurethane foam
26% Cotton
6% Elastane
4% Polyester

Made without natural rubber latex

Pack contents

Compreflex® Standard Calf
1 Compreflex® Standard Calf
1 pair of Liners Basic without compression
1 set Accutab®

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SIGVARIS Britain LTD.
1 Imperial Court
Magdalen Lane
Andover, SP10 5NT

SIGVARIS and SIGVARIS GROUP are registered trademarks of SIGVARIS AG,
CH-9014 St. Gallen/Switzerland, in many countries worldwide.

© 2023 Copyright by SIGVARIS AG, CH-9014 St. Gallen/Switzerland

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EN Instructions for use

Product description

The Compreflex® Standard Calf is an inelastic product that provides precise compression ranges of 21–46 mmHg thanks to the use of the Accutab®. The product is combined with a pair of basic liners to protect the skin.

Intended use

Inelastic wraps are compression therapy devices used on the lower leg. Compression therapy improves venous return, which reduces venous edema. It is primarily used in the initial decongestion phase of pronounced venous edema as an alternative to wrapping with bandages.

This medical device is intended for multiple use and non-sterile.

Effect of medical compression

Compression therapy is used to improve quality of life in chronic venous diseases, for the prevention and initial therapy of venous edema and venous skin changes, and for prevention and therapy of venous leg ulcers. Carefully controlled pressure on the tissue and veins improves the blood return in the legs and increases lymphatic flow. At the same time, it improves microcirculation and therefore the supply to the tissue. Proximal decreasing pressure along the leg must be ensured.

Target group

The intended target group includes patients who have been prescribed a compression product for individual compression therapy by the responsible physician or trained health professional. The target group also includes user groups from health care professions that apply the products to patients.

Indications & contraindications

Indications

Compreflex® Standard Calf

- Edema in chronic venous disease (CVD) due to obstruction or valvular insufficiency
- Chronic venous insufficiency (CVI) including ulcer due to venous stasis

Contraindications

- Severe cardiac insufficiency (NYHA III–IV)
- Phlegmasia coerulea dolens
- Severe PAOD with any of the following: ABI <0.6; ankle pressure <60 mmHg; toe pressure <30 mmHg
- Severe diabetic neuropathy with sensory disorders or microangiopathy with a risk of skin necrosis
- Compression of an existing arterial bypass
- Allergy against compression material

Relative contraindications for which the benefits and risks must be considered, in particular:

- PAOD with ABI between 0.6 and 0.9 and/or ankle pressure between 60 and 90 mmHg
- Patients with compensated heart failure (NYHA I–II) and venous or lymphatic edema
- Diabetes mellitus
- Skin diseases with massive exudation or mixed arterio-venous ulcer
- Compression products to reduce inflammation, pain, or edema in leg erysipelas or cellulitis should only be used in combination with antibacterial treatment.

Risks and side effects

Inelastic wraps can cause skin necrosis and pressure damage to the peripheral nerves if handled improperly. Sensitive skin may develop itching, scaling, and signs of inflammation under the compression devices, so adequate skin care is advised. Regular medical check-ups are required to maintain successful therapy. If acute pain or skin irritation occurs while wearing the compression product, seek medical advice immediately.

Please note: consult your physician. The precautions and contraindications listed here are guidelines only and do not replace medical advice.

Special instructions

- Protect the compression product from ointments, oils, and grease to preserve its service life.
- Compression products must not be worn over open wounds. Wounds should be covered with appropriate dressings.
- Failure to comply with this instruction may reduce the safety and effectiveness of the product.
- All serious incidents that occur when using the product should be reported to the manufacturer and the competent authority.
- Compression is increased when liners with compression pressure are used under a wrap.

Application

Donning the Compreflex® Standard Calf

1. Pull on the provided liner.
2. Before applying, unfold each hook fastener on the product to its full length, roll the ends back, and secure them there so that they cannot get tangled. Slide the foot into the opening between the elastic insert and the product and pull the product up to just below the knee. The textile label is located on the top at the back of the leg and faces outward.