Disposal

There are no special requirements for disposing of the medical compression product. They can be disposed of as residual or household waste. Please comply with local regulations for the disposal of textiles.

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SIGVARIS Britain LTD.
1 Imperial Court
Muggellan Lane
Andover, SP10 5NT

SIGVARIS and SIGVARIS GROUP are registered trademarks of SIGVARIS AG,
CH-9014 St. Gallen/Switzerland, in many countries worldwide.

© 2023 Copyright by SIGVARIS AG, CH-9014 St. Gallen/Switzerland

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EN Instructions for use

Product description

The Compreboot® Standard provides up to 40 mmHg of foot compression and is made of soft, supple neoprene fabric.

Intended use

Inelastic wraps are compression therapy devices used on the foot. Compression therapy improves venous return, which reduces venous edema. It is primarily used in the initial decongestion phase of pronounced venous edema as an alternative to wrapping with bandages. This medical device is intended for multiple use and non-sterile.

Effect of medical compression

Compression therapy is used to improve quality of life in chronic venous diseases, for the prevention and initial therapy of venous edema and venous skin changes, and for prevention and therapy of venous leg ulcers. Carefully controlled pressure on the tissue and veins improves the blood return in the legs and increases lymphatic flow. At the same time, it improves microcirculation and therefore the supply to the tissue. Proximal decreasing pressure along the leg must be ensured.

Target group

The intended target group includes patients who have been prescribed a compression product for individual compression therapy by the responsible physician or trained health professional. The target group also includes user groups from health care professions that apply the products to patients.

Indications & contraindications

Indications

Compreboot® Standard

- Edema in chronic venous disease (CVD) due to obstruction or valvular insufficiency
- Chronic venous insufficiency (CVI) including ulcer due to venous stasis

Contraindications

- Severe cardiac insufficiency (NYHA III–IV)
- Phlegmasia coerulea dolens
- Severe PAOD with any of the following: ABI <0.6; ankle pressure <60 mmHg; toe pressure <30 mmHg
- Severe diabetic neuropathy with sensory disorders or microangiopathy with a risk of skin necrosis
- Compression of an existing arterial bypass
- Allergy against compression material

Relative contraindications for which the benefits and risks must be considered, in particular:

- PAOD with ABI between 0.6 and 0.9 and/or ankle pressure between 60 and 90 mmHg
- Patients with compensated heart failure (NYHA I–II) and venous or lymphatic edema
- Diabetes mellitus

Risks and side effects

Medical compression products can cause skin necrosis and pressure damage to the peripheral nerves if handled improperly. Sensitive skin may develop itching, scaling, and signs of inflammation under the compression devices, so adequate skin care is advised. Regular medical check-ups are required to maintain successful therapy. If acute pain or skin irritation occurs while wearing the compression product, seek medical advice immediately.

Please note: consult your physician. The precautions and contraindications listed here are guidelines only and do not replace medical advice.

Special instructions

- Protect the compression product from ointments, oils, and grease to preserve its service life.
- Compression products must not be worn over open wounds. Wounds should be covered with appropriate dressings.
- Failure to comply with this instruction may reduce the safety and effectiveness of the product.
- All serious incidents that occur when using the product should be reported to the manufacturer and the competent authority.
- Compression is increased when liners with compression pressure are used under a wrap.

Application

Donning the Compreboot® Standard

1. Make sure that the hook fasteners are open.
2. Place your heel in the heel area of the wrap. The textile label attached to the outside in the heel area is used as a guide.
3. Close the hook fasteners around your foot.
4. The hook fasteners may need to be readjusted several times to achieve a comfortable fit.
5. Readjust the product as needed to ensure proper compression.


Doffing the Compreboot® Standard

To remove the product, loosen the hook fasteners starting at the ankle, roll their ends backwards and fasten them there. Once all the hook fasteners have been secured, remove the product from the foot and take off the liner.

Storage, wearing period, and service life

Unless otherwise prescribed by your physician, the compression product can be worn during the day and overnight. However, remove the product once a day to clean the skin and care for any wounds. The product should then be reapplied according to the instructions for use. How often the product needs to be removed for skin checks, hygiene, and/or wound care may vary according to the individual needs of patients and the prescriber's recommendations. A liner should be worn under every compression product. Wearing and washing compression products every day may cause the medical pressure gradient and elasticity of your wrap to diminish. The shelf life date is indicated on the packaging with an hourglass symbol. We guarantee maximum effect for the subsequent six months of wearing.

Care and storage

Follow the instructions for care on the packaging and/or the textile label. Close all hook fasteners before washing. Wash the product in a mesh laundry bag on gentle cycle using a mild detergent. Do not use bleach or fabric softener. Dry the product on low/no heat in a dryer or air dry (air drying is recommended to extend wrap life). Store the product at room temperature and in a dry place away from light. Keep away from direct heat.