# Exadrop®


## Exadrop® SafeSet

---


en	Instructions for use	4
de	Gebrauchsanweisung	5
bg	Инструкция за употреба	7
es	Návod k použití	8
da	Brugsanvisning	10
el	Οδηγίες χρήσης	12
es	Instrucciones de uso	13
et	Kasutusjuhised	15
fi	Käyttöohje	16
fr	Mode d'emploi	18
hr	Upute za uporabu	20
hu	Használati utasítás	21
it	Istruzioni per l'uso	23
kk	Пайдалану нұсқаулығы	25
ko	사용 지침	26
lt	Naudojimo instrukcija	28
lv	Lietošanas instrukcija	29
nl	Gebruikersinformatie	31
no	Bruksanvisning	32
pl	Instrukcja użytkowania	34
pt	Instruções de utilização	35
ro	Instrucțiuni de utilizare	37
ru	Инструкция по применению	39
sk	Návod na použitie	40
sl	Navodila za uporabo	42
sr	Uputstvo za upotrebu	44
sv	Bruksanvisning	45
tr	Kullanım Kılavuzu	47
uk	Інструкції для застосування	48
vi	Hướng dẫn sử dụng	50
zh	使用说明	51

---

**B | BRAUN**

# Exadrop®

## en Instructions for use

### Device description

Single-use, sterile flow control device with / for I.V. administration set of infusion / drugs by gravity – indicated by symbol "G" (Gravity) on primary packaging.

Articles equipped with PVC tube:

ID: 3 mm  
OD: 4.1 mm  
Shore A: 74

Articles equipped with PUR tube (incl. UV-protected tubes):

ID: 3 mm  
OD: 4.1 mm  
Shore A: 87

Exadrop® fulfils the standard EN ISO 8536-13 / ISO 8536-13 for graduated flow regulators which enables the setting or adjustment of a certain flow rate.

Flow rate accuracy:

Exadrop®

Setting	Wheel position	Minimum (ml/h)	Maximum (ml/h)
High	"250"	175	325
Medium	"70"	45,5	95,5
Low	"20"	12	28

Exadrop® SafeSet

Setting	Wheel position	Minimum (ml/h)	Maximum (ml/h)
High	"250"	100	400
Medium	"70"	22	130
Low	"20"	10	28

IV administration sets can be used with devices which are compatible or in conformance with Luer connectors (e.g., ISO 80369-7).

IV administration sets are in accordance with EN ISO 8536-4 / ISO 8536-4.

### Sterility

Sterilized using ethylene oxide. Please refer to the label of primary packaging

Do not use if package is damaged.

Visually inspect the sterile barrier system for breaches and integrity prior to use.

### Intended purpose

Exadrop®:

Precision flow control device with I.V. administration set. Only for gravity infusion.

Exadrop® Inline:

Precision flow control device. Only for gravity infusion.

Exadrop® SafeSet:

Precision flow control device with IV administration set. Only for gravity infusion.

### Indication

IV administration sets are used for infusion therapy e.g.: Fluid volume maintenance, fluid volume replacement, infusion of total parenteral nutrition, infusion of medications (e.g. antibiotics, chemotherapeutic agents) with bottles and bags according to the SPC (Summary of product characteristics) of the drugs / solutions.

### Patient population

The device can be used for all patients for which infusion therapy is prescribed. No gender or age-

related limitations. Exadrop® can be used for adults, pediatric and neonates.

### Intended User

To be used by country-specific authorized persons e.g.:

Healthcare professionals licensed and / or certified physicians, healthcare practitioners, care givers after instruction of a healthcare professional (according to local regulations).

Exadrop® should only be used by healthcare worker who have been adequately trained within the educational background in this technique. Upon medical assessment and adequate instruction, patients and / or caregivers may also be allowed to take over defined handling steps according to national guidelines. Healthcare workers are responsible for ensuring that patients and / or caregivers are instructed in the correct use of Exadrop® and that the instruction of caregiver and / or patient is documented by the healthcare worker according to local documentation guidelines.

### Contra-indications

- IV administration set shall not be used in case of proven incompatibility of the drug / IV solution to one of the used materials. SPCs of the drugs / solutions for contraindications should be considered.
- No administration of blood products or blood components.
- Do not use for pressure infusion > 0,5 bar (50 kPa).

### Side effects

General risk occurring during infusion therapy including air embolism, local- to systemic infection, particulate contamination, medication error including over and under supply. Kinking of the line may result in decrease or stop of infusion fluid. Leakage of drugs / fluids can be an additional hazard to patients, visitors or staff.

### Precautions

- Do not use with pressure infusion devices (infusion pumps).
- Visual check of the device for damages before use.
- Do not use if protective caps are missing or loose.
- Please also consider the corresponding SPCs of the drugs / solutions regarding possible incompatibilities of drugs / solutions.
- Mind the general guidelines for parallel infusions.
- Monitor all connections for tightness over the whole infusion procedure.
- The figures on the scale (= app. ml / h) are only approximate ones. They only refer to the infusion of 0.9 % NaCl solution administered by means of an Vasofix® Braunüle® G 18 being approximately 100 cm static difference of level without CVP. Therefore, the figures on the scale must always be controlled by verification. Otherwise, there is a risk of over- or underdosage for the patient with all subsequent consequences.
- Due to the structure of the fluid filter of Exadrop® SafeSet which achieves a particle retention with a filter efficiency > 80% for a particle size of 3µm, the infusion of highly viscous solutions, lipid emulsions >10% and certain suspensions can affect the flow rate. With regard to this, please also follow the information in the relevant summary of product characteristics.
- Use aseptic techniques in accordance to national guidelines and/or hospital protocols for variants with Safeflow.

### Warning

- Do not re-use. Re-use of single-use devices creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability. Contamination and/or limited functionality of the device may lead to injury, illness or death of the patient.
- Administration of light sensitive drugs requires UV-protected IV administration sets.
- Always check, if incompatibilities with other drugs or with materials is stated in the SPC (Summary of product characteristics) of the drug.
- Do not re-attach the protective cap of the patient connector.
- Do not re-sterilize.
- Do not use after expiration date.

Variants with Safeflow:

- Only use connection devices in compliance with ISO 80369-7. Syringes or connectors which are not in compliance with ISO 80369-7 / EN ISO 80369-7 can damage the swabable valve. Syringes and male luer connectors have a large variety of configurations and can vary significantly in design and dimensions.
- Do not use, in any situation, needles in the needle-free injection port.

### Operating Instructions

1. Open the single package and take out the product.
2. Unfold the line & inspect the IV-Set or the graduated flow regulator in-line set and check whether the slide clamp and flow regulator is opened. Check whether the green air vent cap is closed (only for IV-Set with graduated flow regulator). (Fig. 1 & 2)
3. Connect to accessories (e.g. 3-way stopcocks, manifolds, injection sites) or connect graduated flow regulator in-line set to IV-Set.
4. Remove protective cap from spike (only for IV-Set with graduated flow regulator)
5. Insert the spike into the container (only for IV-Set with graduated flow regulator) Insert administration set with closed snap vent in fluid container. (Fig. 3)
6. Close Flow regulator or graduated flow controller make sure the slide clamp is open. Rotate wheel of graduated flow controller towards "OPEN" in order to release initial position. Set graduated flow controller to "OFF". (Fig. 4)
7. Hang in patient connector in hook of flow regulator housing (only for IV-Set with flow regulator).  
The loop on the Exadrop®-handle is used as tube holding device when the connection of the IV administration set to the IV cannula cannot be established immediately.
8. Hang the container on the stand.
9. Set up fluid level in drip chamber (Fig. 5).
10. Open Flow regulator (if applicable) and graduated flow regulator.  
Set graduated flow regulator to "OPEN". (Fig. 6)
11. Open venting cap (only applies to glass bottles and device variants with venting flaps) (Fig. 7 & 8).
12. Take patient connector from hook and prime the IV line (including accessories (e.g. stopcocks or injection sites) and graduated flow regulator) – air free (Fig. 9).
13. If applicable, follow special priming procedures for accessories e.g. for IV Filter and Flushing sets.

14. Close the flow regulator and / or the graduated flow regulator and venting cap (for glass bottle use only and for device variants with venting flaps), (Fig. 10)
15. Connect to prepared and primed accessories if applicable (e.g. extension line, flushing sets secondary line).
16. Remove patient connector / tubing from flow regulator (only for IV-Set with flow regulator).
17. Remove protective cap from patient connector.
18. Lock patient connector to the venous access (Fig. 11).
19. Open Flow Regulator (if applicable) and/or graduated flow regulator and venting cap (for glass bottle use only) and start infusion (Fig. 12).
20. Adjust flow rate by turning the wheel and monitor infusion rate.  
Set desired flow rate on the markings of the wheel by working back from "OPEN".
21. Control the drop rate during infusion.
22. Important note:  
When the prepared IV-system is laid down horizontally, e.g. on a tablet, or the spiked glass container is placed on a table, it is necessary to close the air vent cap beforehand! Only re-open the air vent cap when the infusion is running.  
Always ensure to establish a 20 cm siphon (Fig. 12 Et 13).
23. Stop of therapy: interrupt infusion by closing flow regulator or by closing slide clamp or by closing graduated flow regulator.  
  
With the slide clamp above the control the infusion can be interrupted without setting back the control device to "OFF". After having opened the clamp the flow rate remains the same if the pressure conditions in the infusion system have not changed in the meantime. (Fig. 14)

#### Duration of use

Duration of use is dependent on the intended therapy according to the SPC of the drug or solution. In general, change of administration sets should occur according to national guidelines (e.g. CDC guideline) and / or hospital protocols.

#### Disposal

Disposal acc. to local guidelines and / or clinical protocols.

If applicable use integrated spike protection at flow regulator housing for disposal.

#### Storage and handling conditions

Storage conditions:

Keep away from sunlight.

Keep dry.

#### Notice to the user

If, during the use of this product or as a result of its use, a serious incident has occurred, please report it to the manufacturer and / or its authorized representative and to your national authority.

If further Instructions for Use are required, they can be requested from the manufacturer or can be obtained from the B. Braun homepage: <https://eifu.bbraun.com/>.

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B. Braun Melsungen AG