
klinidur® cohesive bandage

INSTRUCTIONS FOR USE

KLINION® WOUND CARE

IFU_KWC062_01
20-06-2023 / 0_030102

ENGLISH

Product Description

Klinidur Cohesive bandage is an elastic, beige colored, cohesive short stretch compression bandage.

Made of 100% cotton with a latex-free cohesive coating on both sides.

Not made with natural rubber latex. Medical device.

Indications

Klinidur Cohesive bandage is indicated for the treatment of venous leg ulcers from the lower extremity, lymphatic oedemas and other clinical conditions where compression is appropriate.

Contra-indications

- Do not apply directly on open wounds or breached skin.
- Decompensated heart insufficiency
- Septic phlebitis
- Advanced diabetic microangiopathy
- Phlegmasia cerulea dolens
- Severe arterial occlusive disease in patients with an ankle-brachial pressure index (ABPI) of less than 0.8
- Known hypersensitivity to any of the materials of the components.

Directions for use

Application

1. Apply the bandage with the foot in dorsiflexion (90°). Ensure that this angle is maintained through the entire application.
2. Start by applying the bandage to the metatarsophalangeal joint in a circular manner, wrapping two turns to –fix the bandage. Use moderate tension and slightly elevate the lateral edge of the foot.
3. Wrap the bandage over the back of the foot and secure the heel by using a figure of eight turns with minimum stretch. Ensure full coverage of the heel.
4. Continue wrapping the bandage in a circular manner up the leg with a minimum stretch and an overlap of 50%. Finish bandaging around 2 cm below the popliteal space. Cut off the excess bandage, and if necessary secure with adhesive tape. Carefully press down the bandage to ensure full cohesion.

Removal and disposal

Klinidur Cohesive bandage can be removed with bandage scissors or by unwrapping.

Follow your current national legislation, standards and guidelines, which regulate the disposal of medical waste for the disposal of the compression system.

Precautions

- For single use only.
- Do not reuse. Reuse can compromise product integrity or lead to device failure.
- Always apply under the supervision of a healthcare professional. It is recommended to have an adequate training for the correct application of the bandage.
- Check that the patient's ABPI is higher than 0,8 before applying the bandage.
- Patients with known arterial insufficiency, decompensated heart insufficiency, or advanced diabetic microangiopathy may not tolerate compression; it should only be used under medical prescription, with careful monitoring and after appropriate treatment.
- Wrapping too tightly may impair circulation. Monitor the area of application frequently for signs of discoloration, pain, numbness, tingling or other changes in sensation and swelling. If these symptoms appear, patients should remove the bandage immediately and inform the healthcare professional.
- Before applying, it is necessary to perform the usual diagnostic examinations. Wounds must be cleansed and covered with wound dressings appropriate to the wound conditions.
- Follow current international guidelines and local regulations for infections control.
- If the bandage becomes wet, loose fitting or when it no longer conforms to the shape of the bandaged area, it should be replaced.

Storage

The products should be stored under ambient conditions (5°C-30°C), in the original packaging. In a cool, dry and well ventilated place, away from direct sunlight.

In case a serious incident has occurred in relation to the device: Inform Farmaban S.A. as listed on the packaging, your local distributor and your competent healthcare authority.

 **Manufactured by**
Farmaban S.A.
Carrer Vil·la de Claret, 6
08272 Sant Fruitós de Bages
Barcelona, Spain
www.farmaban-sa.com

 **Distributed by**
Medeco B.V.
Brandpuntlaan Zuid 14
2665 NZ Bleiswijk
The Netherlands

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