Baxter

Hemopatch


SEALING HEMOSTAT

Resorbable Collagen-Based Sealing Hemostat

REF 1506257RT (5 units)

REF 1506256RT (3 units)

REF 1506253RT (3 units)

Contents

(en)	HEMOPATCH Sealing Hemostat	1	(it)	HEMOPATCH Sigillante emostatico	53
(bg)	Запечатващ хемостатик HEMOPATCH	5	(lt)	HEMOPATCH hermetizuojanti hemostazinė priemonė	57
(cs)	Těsnící hemostat HEMOPATCH	9	(lv)	HEMOPATCH hermetizējošs hemostatisks līdzeklis	61
(da)	HEMOPATCH forseglende hæmostatikum	12	(nl)	HEMOPATCH Sealing hemostaticum	65
(de)	HEMOPATCH Hämostatikum mit versiegelnden Eigenschaften	16	(no)	HEMOPATCH forseglende hemostatikum	69
(el)	Αιμοστατικό στεγανοποίησης HEMOPATCH	21	(pl)	Hemostatyk uszczelniający HEMOPATCH	73
(es)	Sellante hemostático HEMOPATCH	25	(pt)	Selante hemostático HEMOPATCH	78
(et)	HEMOPATCH'i hermetiseeriv hemostaatikum	29	(ro)	HEMOPATCH plasture hemostatic	82
(fi)	HEMOPATCH tiivistävä hemostaatti	33	(ru)	Герметизирующее гемостатическое изделие HEMOPATCH	86
(fr)	Patch hémostatique d'étanchéité HEMOPATCH	37	(sk)	Hemostatikum na zastavenie krvácania HEMOPATCH	90
(hr)	HEMOPATCH hemostatičko ljepilo	41	(sl)	Lepilni hemostatik HEMOPATCH	94
(hu)	HEMOPATCH lezáró haemostaticum	45	(sv)	HEMOPATCH förseglande hemostas	98
(is)	HEMOPATCH blæðingahemjandi sárarúði	49	(tr)	HEMOPATCH Sızdırmaz Hemostat	102

HEMOPATCH Sealing Hemostat Resorbable Collagen-Based Sealing Hemostat

DESCRIPTION

HEMOPATCH Sealing Hemostat ("HEMOPATCH") consists of a soft, thin, pliable, flexible pad of collagen derived from bovine dermis, coated with NHS-PEG (pentaerythritol polyethylene glycol ether tetra-succinimidyl glutarate).

Due to its flexible structure, the application of HEMOPATCH to the site where hemostasis / sealing is desired is easily controlled. For differentiation, the non-coated side is marked with blue squares using a biocompatible colorant.

CLINICAL BENEFITS

- Seals tissue, stops and prevents leaks
- Fast effective hemostasis within 2 minutes
- Safe and effective sealing hemostat
- Ready-to-use (no preparation required)
- Resorbed and replaced by host tissue

HEMOPATCH is supplied in the following 3 sizes:

- 27 x 27 mm **REF** 1506257RT (5 units)
- 45 x 45 mm **REF** 1506256RT (3 units)
- 45 x 90 mm **REF** 1506253RT (3 units)

INTENDED PURPOSE

HEMOPATCH is an absorbable collagen pad intended for sealing and hemostasis.

INDICATION

HEMOPATCH is indicated as a hemostatic device and surgical sealant across various soft tissues (cardiovascular, connective tissue, parenchyma, serosa, viscera), and dura for procedures in which control of mild or moderate bleeding or leakage of other body fluids or air by conventional surgical techniques is either ineffective or impractical. HEMOPATCH may be used to augment dura closure techniques to close small dural defects ($\leq 3\text{mm}$) following traumatic injury, excision, retraction or shrinkage of the dura mater.

INTENDED USER

HEMOPATCH is intended to be used by licensed healthcare professionals/surgeons and clinicians qualified by proper education in a medical setting.

INTENDED PATIENT POPULATION

Adult and pediatric patients of any gender undergoing a surgical procedure. Only limited clinical data is available for the pediatric patient population.

MECHANISM OF ACTION AND PERFORMANCE CHARACTERISTICS

The performance of HEMOPATCH involves the interaction of two components that are able to create rapid and lasting hemostasis by sealing off the bleeding surface.

The white, tissue-facing side of HEMOPATCH is covered with a thin layer of NHS-PEG. The NHS-PEG, when in contact with blood or other body fluids forms a hydrogel which enhances its adhering properties and seals the tissue surface.

In contact with blood, collagen induces aggregation of platelets. Platelets deposit in large numbers on the collagen structure, degranulate, and release coagulation factors that, together with plasma factors, enable the formation of fibrin. The structure of HEMOPATCH provides a three dimensional matrix that readily absorbs fluids, and provides additional mechanical strengthening of the clot. Preclinical studies suggest that resorption and replacement with host tissue occurs in approximately 6 – 8 weeks with little tissue reaction. The rate of resorption and replacement with host tissue can be longer as it is dependent on several factors, including the amount of product left in situ and the site of use.

CONTRAINDICATIONS

Do not compress HEMOPATCH into blood vessels or use intravascularly. The device must not be used in patients with known hypersensitivity to bovine proteins or brilliant blue (FD&C Blue No. 1 [Blue 1]).

WARNINGS

HEMOPATCH is not recommended to be used in pulsatile, severe bleedings. The use of HEMOPATCH is not recommended in the presence of an active infection.

When used in, around, or in proximity to foramina in bone, areas of bony confine, the spinal cord, the brain and/or cranial nerves, care should be exercised to avoid overpacking as HEMOPATCH may expand upon absorption of liquid, creating the potential for neural damage.

HEMOPATCH is not intended as a substitute for meticulous surgical technique and the proper application of ligatures or other conventional procedures for hemostasis and sealing.

PRECAUTIONS

For single use only. Do not re-sterilize. Reuse of single-use devices creates a potential risk of infection.

Do not store without aluminum blister.

Discard HEMOPATCH if the packaging is damaged.

Do not apply the surface marked with blue squares to the application site.

Do not wet HEMOPATCH prior to application.

HEMOPATCH is self-adherent and is not intended to be sutured.

Do not apply HEMOPATCH on a dry tissue surface or lesion, or in the presence of any excessive collection of body fluids or irrigation solution. NHS-PEG only forms an adhering hydrogel on a moist surface, that is neither too dry nor covered with excessive fluid collections. HEMOPATCH must be in direct contact with the wound or tissue lesion moistened with body fluid such as blood or lymphatic, or with sodium bicarbonate solution¹ (concentration between 4.2% to 8.4%).

When HEMOPATCH has adhered to tissue, it should only be repositioned per the surgeon's discretion, as this may result in damage to the tissue.

Do not use HEMOPATCH on bone surfaces where adhesives, such as methylmethacrylate or other acrylic adhesives, will be required to attach a prosthetic device. Microfibrillar collagen has been reported to reduce the strength of methyl methacrylate adhesives used to attach prosthetic devices to bone surfaces.

Applying fibrin or other sealants to target tissue prior to application of HEMOPATCH could result in suboptimal performance of the product.

The safety and performance of HEMOPATCH has not been established in pregnant and/or lactating women.

ADVERSE EVENTS

The following adverse events have been identified during the use of HEMOPATCH: Seroma formation, Hematoma, Potentiation of infection, Inflammatory reaction, Leakage.

These reactions are reported voluntarily from a population of uncertain size and it is not always possible to reliably estimate their frequency or establish a causal relationship to device exposure.

Collagen-Based Agents Adverse Events:

Adverse events which may be possibly or probably related to the use of any collagens include: Seroma formation, Hematoma, Potentiation of infection, Inflammatory reaction, Foreign body reaction, Adhesion formation, Allergic reaction, Nerve compression, Leakage.

Serious Adverse Events Reporting

For a user and/or patient if, during the use of this device, or as a result of its use, a serious incident has occurred, please report this incident to the manufacturer, and/or its authorized representative, and to the competent authority of the Member State in which the user and/or patient is established.

DRUG INTERACTIONS

Chlorine-releasing antiseptic agents as well as tannins and caustics should not be used in combination with collagen at the application site.

PREPARATION

HEMOPATCH comes ready to use in sterile packages and must be handled accordingly. Use only undamaged packages. Inspect the sterile barrier (aluminum blister) integrity prior to use. Single use only. Do not re-sterilize. HEMOPATCH is MR safe.

Circulating Nurse:

- Open outer, non-sterile package (aluminum blister) and,
- Present the inner, sterile blister pack to the scrub nurse.

Scrub Nurse:

- Deliver the inner, sterile package into the sterile field, using aseptic technique. Once removed from the outer blister HEMOPATCH can be used throughout surgery provided it remains dry and sterile.
- Open the inner, sterile package and present the device to the surgeon.
- HEMOPATCH should be maintained dry prior to application. Do not soak in solution prior to application (e.g. physiological saline solution, antiseptic agents).

¹ Footnote: ¹the grade of sodium bicarbonate solution to be suitable for parenteral human use.

- Dry gloves and instruments (e.g. forceps, scissors) should be used to handle HEMOPATCH.
- In situations in which HEMOPATCH has been removed for application from the inner blister, but application is delayed, it can be used throughout surgery, provided it remains dry and sterile.

ADMINISTRATION

- **Do not apply intravascularly.**
- Select the appropriate size of the pad so that it overlaps the margins of the bleeding surface or tissue lesion by about 1 cm. The pad may be cut to the desired size and shape.
- Dry gloves and surgical instruments (forceps, scissors) should be used to handle, cut and apply HEMOPATCH.
- Multiple pads may be used for larger bleeding surfaces or tissue lesions. When using multiple pads, ensure a 1 cm overlap where the pads contact each other.

Maximum number of patches to be used

HEMOPATCH size	Adults	Children (≥ 10 Kg)
45 x 90 mm	7	1
45 x 45 mm	14	2
27 x 27 mm	28	4

METHOD OF APPLICATION

- Before application of HEMOPATCH, aspirate any excessive collection of body fluids such as blood or lymphatic or irrigation solutions to avoid floating of HEMOPATCH and ensure immediate direct contact with target tissue.
- Apply the dry HEMOPATCH with the non-marked white surface in direct contact with the wound or tissue lesion moistened with body fluid. In the absence of such wound fluids, sodium bicarbonate solution² (concentration between 4.2% to 8.4%) can be used to moisten the tissue prior to application of HEMOPATCH.
- After application of HEMOPATCH, immediately use a dry gauze or patty and hold in place with gentle, uniform pressure over the entire

pad surface for 2 minutes. When using HEMOPATCH laparoscopically, after trocar passage, dry graspers (with gauze, if feasible) can be used to place and approximate HEMOPATCH with uniform pressure for 2 minutes directly on the target site.

- When applying HEMOPATCH, minimize contact with bloody or moist surgical instruments, gauzes or gloves due to the affinity of collagen to blood and body fluids.
- Gently remove the gauze or patty from the pad. Gentle irrigation may also help in removing the gauze or patty without dislodging HEMOPATCH from the surgical site.
- Leave HEMOPATCH in situ after hemostasis or sealing has been achieved.
- If hemostasis or sealing has not been satisfactorily achieved, HEMOPATCH may be removed at the surgeon's discretion. A new pad can be applied according to the steps above.
- When HEMOPATCH has adhered to tissue, it should only be repositioned per the surgeon's discretion, as this may result in damage to the tissue.

SHELF LIFE AND STORAGE

Store in a dry place at 2° to 30°C.

Expiry date: See imprint on outer package.

Do not use after the expiry date printed on the label.

INSTRUCTIONS FOR COMPLETION OF IMPLANT CARD*

Enter the following information into the spaces indicated on the Implant Card

- Name and address of the Healthcare Institution or provider

- Name of patient or patient Identification

- Date of implantation

² Footnote: ²the grade of sodium bicarbonate solution to be suitable for parenteral human use.

DISPOSAL

Any packaging, or a HEMOPATCH exposed to human tissue/fluids should be disposed of per hospital protocols. Any open, unused HEMOPATCH should be discarded per hospital protocols if sterility is compromised.


ADDITIONAL INFORMATION FOR PATIENT AND HEALTHCARE PROFESSIONALS

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE (SSCP)*


For pertinent safety and clinical information about this product, the Summary of Safety and Clinical Performance intended for healthcare professionals (Part I) and patients (Part II) can be found at the following location: <https://ec.europa.eu/tools/eudamed>

Until SSCPS are available on the EUDAMED website, you can find the SSCP at <http://sscp.baxter.com>

*for patients residing in the European Union

Definition of Symbols			
	Medical Device		Temperature limit
	Consult instructions for use		Use-by date
	Caution		Date of manufacture
	Do not re-use		Batch code
	Sterilized using irradiation		Catalogue number
	Do not re-sterilize		Unique Device Identifier

Definition of Symbols

	Single sterile barrier system with protective packaging inside		Manufacturer
	Do not use if package is damaged <i>and consult instructions for use</i>		Importer
	Contains biological material of animal origin		Authorized Representative in the European Community
	Not made with natural rubber latex		Patient information website

The electronic version of the Instructions for Use can be found at the following location: <http://edocs.baxter.com>

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Baxter and Hemopatch are trademarks of Baxter International Inc. or its subsidiaries.
Baxter Deutschland GmbH