EN | Hyaluronic acid sodium salt solution for intra-articular use vials and pre-filled syringes

Product information

REPLASYN is composed by a viscous solution of hyaluronic acid sodium salt (500-730 kDa ca.), in buffered physiological sodium chloride, pH 6.8-7.5. The solution (20mg/2ml) is contained in vials/syringes, is provided sterile and it is free from bacterial endotoxins.

Hyaluronic acid, a polysaccharide of the glycosaminoglycan family, is naturally present in many human tissues such as cartilage and synovial fluid. Endogenous hyaluronic acid provides viscosity and elasticity to synovial fluid, which are fundamental properties for joint function. In traumatic and degenerative joint disorders, an insufficient amount of hyaluronic acid and a lower degree of polymerization occur in synovial fluid and cartilage resulting in an impairment of joint function and in a painful symptomatology. Extensive data from literature indicate that intra-articular administrations of exogenous hyaluronic acid sodium salt are capable to restore the rheological properties of synovial fluid. REPLASYN is a treatment for temporary synovial fluid replacement in patients affected by degenerative or traumatic osteoarthritis, which causes an alteration of the functional performances of the synovial liquid. Intra-articular injection of REPLASYN at weekly intervals, improves the symptoms reducing pain (at rest and on motion) and ameliorating joint function.

Composition

Principal component: Hyaluronic acid sodium salt 20mg/2ml
Other components: Sodium chloride, Disodium hydrogen phosphate dodecahydrate, Sodium dihydrogen phosphate dihydrate, Water for injection.

Indications

Temporary synovial fluid supplementation in traumatic and degenerative joint disorders such as osteoarthritis affecting ankle, trapeziometacarpal joint, temporomandibular joint, knee, hip and shoulder. The product is indicated in the treatment of persistent pain after initial failure of analgesics or in case of failure or intolerance to non-steroidal anti-inflammatory drugs.

Dosage

A treatment cycle for osteoarthritic disorders in ankle, trapeziometacarpal joint and temporomandibular joint consists in 5 injections administered at weekly interval.

A treatment cycle for osteoarthritic disorders in knee, hip and shoulder consists in 3 injections administered at weekly interval.

Administration

- The product is intended only for intra-articular injections.
- Product administration should be performed exclusively by qualified physicians (e.g. orthopaedic surgeon, rheumatologist, physiatrist, radiologist, anaesthetist, sports doctors, etc.)”
- All the rules regarding the asepsis and the injection technique should be followed.
- Remove joint effusion, if present, before injecting REPLASYN. Patients should be carefully examined prior to administration to determine signs of acute inflammation, and the physician should evaluate whether REPLASYN treatment should be initiated when objective signs of inflammation are present.
- Disinfectants containing quaternary ammonium salts should not be concomitantly used because hyaluronic acid can precipitate in their presence.
- As with any invasive joint procedure, it is recommended that care should be taken not to overburden the joint immediately after the intra-articular injection.

Contraindications

REPLASYN should not be administered to patients with:

- known hypersensitivity to any component of the product.
- severe liver pathologies.

Intra-articular injections of REPLASYN are contraindicated in cases of infections or skin diseases in the area of the injection site.

Use in specific populations

REPLASYN has not been tested in pregnant or breastfeeding women nor in children under 18 years of age and therefore its use is contraindicated in these patient populations.

Warnings and precautions

- Do not use the product after the expiry date reported on the package.
- Do not use the product once the packaging is open or it is damaged: the sterility of the solution is guaranteed as long as the primary package is closed and undamaged.
- Caution should be used in patients on treatment with anticoagulants and/or antiplatelet agents since the injection procedure could cause haematomas at the site of injection.
- The product is for single use, that means it is intended to be used once only for a single patient. The assembled syringe must be discarded immediately after use, regardless of whether or not the solution has been completely administered.
- If this product is reprocessed and/or reused, Fidia farmaceutici cannot guarantee performance, functionality, material structure, or cleanliness or sterility of the product. Reuse could lead to illness, infection and/or serious injury to the patient or user.
- After use, dispose according to applicable national practice.
- Keep out of reach of children.

Undesirable effects

Arthralgia and local reactions such as pain, swelling/effusion and warmth or redness at the injection site have been reported. Usually such symptoms are transient and disappear spontaneously within a few days by resting the affected joint and applying ice locally. Only sporadically, these events have been more severe and longer lasting. Very rare cases of fever, sometimes associated with local reactions, have been also reported. In some cases, these undesirable effects could be related to pre-existing medical conditions of the patient. As for any intra-articular treatment, septic arthritis may rarely occur when general precautions for injections are not observed or the site of injection is not aseptic. Local or systemic allergic reactions (eg.: rash, pruritus) due to individual hypersensitivity have been sporadically recorded. Other very rare undesirable effects, such as headache, dizziness, nausea and diarrhea, have also been reported.

Storage

The product must be stored at room temperature not exceeding 25°C in the original package, protected from light. Do not freeze.

How supplied

Box containing 1 vial sealed in blister.

Box containing 1 pre-filled syringe sealed in blister.

Box containing 5 pre-filled syringes, each sealed in blister.

The solution of 20 mg/2 ml hyaluronic acid sodium salt contained in the vial/syringe is sterilised using steam.

Date of the latest revision of the instructions for use: December 2020

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