
Suprasorb® H

Hydrokolloid-Verband
Pansement hydrocoloïde
Hydrocolloid Dressing
Apósito de hidrocoloide
Penso hidrocolóide
Medicazione idrocolloidale
Hydrocolloïdverband
Hydrokolloid-forbinding
Hydrokolloidförband
Hydrokolloid-sårforbinding
Hydrokolloidisiide
Hydrokoloidní obvaz
Hydrokoloidný obváz
Hidrokolloid kötszer
Opatrunek hydrokoloidowy
Υδροκολλοειδές επίθεμα
ضمادة غرواني مائي

STERILE R

34°F
1°C

86°F
30°C

T&L Co., Ltd., 70-17, Wonam-ro, Won-gok-myeon, Anseong-si, Gyeonggi-do, 17554, Korea

EC REP KTR Europe GmbH, Mergenthalerallee 77, 65760 Eschborn, Germany

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en Instructions for Use

Product Description

Suprasorb® H Hydrocolloid Dressing has a two layer construction composed of an absorbent hydrocolloid matrix laminated to a breathable polyurethane film backing. Hydrocolloid matrix rapidly absorbs exudate to form a cohesive gel. In addition to excellent absorbency, the breathable polyurethane film backing provides a consistently high rate of moisture vapor transmission. Together, these features ensure an optimal moist wound environment, minimize the chance for damage to healthy peri-wound skin and provide wear time for up to seven days.

Indications

Suprasorb® H STANDARD is primarily indicated for the management of lightly to moderately exuding pressure sores and leg ulcers. Suprasorb® H THIN and Suprasorb® H STANDARD may be used for superficial burns, superficial partial-thickness burns, donor sites, diabetic ulcers, post-operative wounds and skin abrasions.

Instructions for Use

- Clean the wound area with normal saline or water.
- Gently dry the surrounding skin and ensure that it is free of any residue.
- Choose the size and shape of dressing which is most suitable. It should extend 1.5 to 2.0 cm beyond the wound margin.
- Remove the large part of the protective paper from the dressing.
- Apply the exposed part of the dressing, centering it over the wound bed.
- Avoid touching the adhesive surface to minimize the risk of contamination. Then remove the small part of the protective paper from the dressing and apply it.
- Finalize application of the dressing, making absolutely sure that the edges of the dressing adhere well. This initial adhesion will improve as the dressing warms up. If necessary, use medical tape around the edges to avoid roll-up.

Frequency of Change

The frequency of change will vary from the condition of the wound and the amount of the exudates produced. The dressing should be left in place for the longest possible time. In case of leakage the dressing should be changed immediately. Change the dressing after a maximum of 7 days or when full absorbency has been reached.

Removal Procedure

Press down on the skin and carefully lift the dressing starting at one edge. Pull towards the center of the wound until the dressing is removed. Clean the wound if necessary.

Storage Conditions

Store at room temperature, avoiding high humidity and refrigeration.

Precautions

Hydrocolloid wound dressings should only be used under the advice of a health care professional. Inappropriate use, particularly on patients with fragile skin, can result in skin stripping. A medical practitioner should treat wounds with signs of clinical infection before treatment with Suprasorb® H. Ulcers predominantly or entirely caused by arterial insufficiency should be inspected regularly by a medical practitioner.

Contraindications

- Third-degree burns
- Clinically infected wounds
- Wounds at high risk of infection
- Patients with a known allergy to hydrocolloid dressings

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