Suprasorb® P SensiFlex border, multisite border

Silikon-Schaumverband
Pansement hydrocollulaire silicone
Silicone Foam Dressing
Apósito de espuma y silicona
Ligadura esponjosa de silicone
Medicazione in schiuma di silicone
Siliconen schuimverband
Silikone-skumforbinding
Silikonskumförband
Silikon-skumbandasje
Silikonivaahattosidos
Silikonový pěnový obväz
Silikónový penový obväz
Szilikon habkötés
Opatrunek z pianki w warstwą silikonu
Obloga iz silikonске пене
Povoj sa silikonomk pljenom
Απορροές ολικονοφού επίθεμα
Silikonlu köpük yara örtüsü

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en Instructions for use

Product description

Suprasorb® P SensiFlex (multisite) border is composed of:

- a silicone wound contact layer,
 - a flexible wound pad consisting of three layers: an absorbent polyurethane foam layer, a non-woven distribution layer and a super-absorbent core; the distribution layer and the absorbent core are perforated, which makes the wound dressing more flexible,
 - a vapour-permeable and waterproof protective backing.
- The five-layer construction facilitates dynamic fluid management to provide and maintain optimal moist wound environment which leads to the promotion of faster wound closure and may help to reduce the risk of maceration. The gentle silicone layer can be lifted and repositioned without losing its adherence. Also, this layer can minimise pain to patients and trauma to wounds and the surrounding skin during dressing changes. The absorbent padding layer provides superior absorption and locks exudates away.

Product composition

Silicone, polyurethane, polyacrylate, fluff pulp and polyolefin

Intended purpose

The Suprasorb® P SensiFlex silicone foam dressing can be used to treat all wounds, provided they are

 moderately to heavily exuding and

 superficial to

 deep* wounds.

Indications

- Pressure ulcers
- Leg and foot ulcers
- Traumatic wounds
- Surgical wounds
- Skin tears

The dressing may also be used as part of a prophylactic therapy to help prevent pressure ulcers.

Contraindications

Known allergy and/or hypersensitivity to any of the product components

Application instructions

May be applied by lay persons under supervision of health care professionals.

The product can be used under compression bandaging.

1. Prepare the wound: Clean the wound in accordance with current clinical standards. Dry the skin surrounding the wound thoroughly.
2. Apply the wound dressing: a. The size of the dressing depends on the wound size. The wound pad should cover the dry surrounding skin by at least 1-2 cm.
- b. Remove the first release liner and apply the adherent side to the wound.
- c. Remove the remaining release liner(s) and smooth down the border on the skin. Do not stretch the dressing.
- *Where necessary, for deep wounds apply a packing rope as primary dressing (e.g. Suprasorb® A, Suprasorb® A + Ag, Suprasorb® Liquecel).
3. Dressing change: The dressing can remain on the wound for several days, dependent on wound condition and the level of exudate.

Precautions and warnings

- The product is for single use only.
- Do not use it if the packaging is damaged or opened.
- Do not use Suprasorb® P SensiFlex silicone foam dressing together with oxidising agents such as hypochlorite solutions or hydrogen peroxide.
- If reddening or sensitisation occur discontinue use and consult a health care professional for appropriate treatment.
- Because of the lack of data, cautious use is advised in sensitive population groups such as infants, children, pregnant or nursing women.

Note

In case of serious incident, contact the manufacturer and the responsible health authorities.

Storage and disposal

1. Please store the product in a clean and dry place in specified temperature range and away from sunlight.
2. Disposal should be handled according to local environmental procedures.

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