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### Actinol® Solution Plus

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### INSTRUCTIONS FOR USE

#### Actinol® Solution Plus

#### Cleaning, irrigation, moistening and care solution for wounds

Actinol® Solution Plus is a cleaning, irrigation, moistening and care solution used together with materials such as bandages, dressings, pads and wound fillers for acute, chronic and infected wounds as well as 1<sup>st</sup> and 2<sup>nd</sup> degree burns.

Actinol® Solution Plus cleans wounds mechanically. Surface active agents in its composition provide a convenient and painless cleaning needed before wound care by decreasing the surface tension. Actinol® Solution Plus contributes to the wound healing process by ensuring a highly positive cleansing and moisturizing effect on applied areas.

Actinol® Solution Plus can also be used for antiseptic cleansing and softening of dressing materials. Bandages and dressings dried and adhered to the wound can be wetted by Actinol® Solution Plus before removal, so the microorganisms settled on them can be eliminated and they can be easily removed without causing any secondary trauma to the wound. The antimicrobial effect of Actinol® Solution Plus offers an ideal solution to prevent biofilm formation on bandages and dressings.

Actinol® Solution Plus does not inhibit granulation or epithelialization. It has neither irritating nor pain giving properties. It causes no sensitization and can be used safely for a long time. The ingredient polyhexanide is millionfold clinically proven, dermatologically harmless and is considered to be the first choice as wound treatment agent<sup>1</sup>.

<sup>1</sup>For a healthy wound recovery process, wound cleaning and care should be performed in regular intervals.

<sup>1</sup>Consensus recommendation of the German Society for Wound Healing and Wound Management e.V. (DgWfW) 2004 & 2009.

#### Indications

- Infected or uninfected acute wounds (e.g. traumatic wounds that occur as the result of cuts, tears, excoriations, crushing, being bitten, being stung, etc.)
- Infected or uninfected chronic wounds (e.g. diabetic ulcers, venous ulcers, decubitus (pressure sores), etc.),
- Infected or uninfected thermal or chemical wounds (e.g. 1<sup>st</sup> and 2<sup>nd</sup> degree burns, chemical burns, etc.),
- Wounds that caused by electrical current,
- Radiation wounds,
- Entry regions of urological catheters,
- Maintenance of peristomal skin, fistulas and abscesses,
- Suitable for cleaning and maintenance of wound surface prior to dressing.

Actinol® Solution Plus should be used for wound surface cleaning and irrigation before any treatment. Actinol® Solution Plus can be applied in combination with Actinol® Gel Plus. The wound should be cleaned with Actinol® Gel Plus to remove any contaminants and deposits. Bandages, dressings, pads etc. soaked in Actinol® Solution Plus can be used to clean wounds.

Actinol® Solution Plus can also be used for antiseptic cleansing and softening of encrusted dressing materials that are dried and adhered to the wound surface and difficult to remove. Actinol® Solution Plus softens the encrusted bandage and dressing, while eliminating the pathogenic microorganisms that have settled on them and prevents the wound from being infected.

Actinol® Solution Plus can be applied directly onto the wound as well as onto the bandage and gauzes. Materials

## such as bandage, gauze, wiping clothes, etc. should not be used again. Actinol® Solution Plus can be heated to body temperature prior to use.

### Usage with dressings and bandages

Thoroughly wet dressings and bandages adhered to the wound surface with Actinol® Solution Plus, wait ca. 5-10 min. For softening and remove them carefully from wound. In order to maintain a moist wound climate, the dressing can be changed daily, if necessary, several times a day. Actinol® Solution Plus is compatible with all modern wound dressings.

### Catheter care

Possible infections can be avoided by regular cleaning catheter inlet ports, flexible connectors and its environment by Actinol® Solution Plus. Spray Actinol® Solution Plus on the area from a short distance and wait ca. 1-5 min. The ingredients of Actinol® Solution Plus are particularly mucous membrane-friendly, cleanse catheter and PEG material without damaging it. As with wound care, encrusted sticking plasters can be removed easily after wetting with Actinol® Solution Plus.

Side Effects Actinol® Solution Plus contains polyhexanide. Even though it is not irritating, it can cause allergic reactions in the bands and the patient's hands. Actinol® Solution Plus is reported that polyhexanide may cause an anaphylactic shock in rare cases (less than 1 out of 10,000). When allergic skin testing is recommended should not use. Patients hypersensitive to polyhexanide should not use the product. The biocompatibility index of Actinol® Solution Plus is found to be >1 and has a high tissue compatibility. Actinol® Solution Plus may cause allergic reactions such as itching (urticaria) and rash (eczema). In very rare cases, a mild burning sensation may occur after application of Actinol® Solution Plus, however this sensation disappears in a couple of minutes.

**Contraindications** Actinol® Solution Plus in the following cases: It is not suitable for use in the cartilage and aseptic joint surgeries, in the middle and internal ear and in the eye (Actinol® Solution Plus contacts the aseptic cartilage, wash immediately with Ringer solution or physiological saline solution).

- If the patient is known or suspected to be hypersensitive to any of the ingredients,
- On central nervous system or meningis,
- Along with anionic surfactants, cleansing soaps, creams, oils, enzymes, etc. These substances should be removed fully from the wound surface before use.
- For peritoneal lavage and washing,
- It should not be used in the first four months of pregnancy. Depending on the risk-benefit ratio, it should be used later in pregnancy should be evaluated according to the doctor's recommendation.

**Warnings and Precautions** Read the instructions carefully before use. Only for external application onto skin. Do not use for infusion or injection purposes. Do not swallow! Keep out of reach of children. Do not use if the sealing under the cap is torn or broken. Close the bottle immediately after use to avoid any contamination. The rim of the bottle should be protected from direct contact to avoid contamination while using. Bottles, which directly contact the wound or are contaminated in some other ways should be disposed. Keep only in original packaging contact from sunlight. Dispose according to local regulations.

**For newborns and infants:** Due to the insufficiency of clinical information, Actinol® Solution Plus should be used only under a strict and close medical supervision on newborns and infants.

This product is a Class II Medical Device according to the 93/42/EEC and the 2007/47/EC amendments.

**Ingredients:** Poloxamer, Polyhexamethylene biguanide (PHMB, 0.1%) and excipients.

**Appearance and odor:** Clear, colorless and almost odorless solution.

**Shelf life:** Please see the expiry date. Close the cap tightly and store in a cool, dry and frost-free environment. Shelf life is 2 (two) years at room temperature. Do not store an opened bottle above 35 °C. Consume the open product within 8 (eight) weeks.

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Actinol® Solution Plus should be used for wound surface cleaning and irrigation before any treatment. Actinol® Solution Plus can be applied in combination with Actinol® Gel Plus. The wound should be cleaned with Actinol® Gel Plus to remove any contaminants and deposits. Bandages, dressings, pads etc. soaked in Actinol® Solution Plus can be used to clean wounds.

Actinol® Solution Plus can also be used for antiseptic cleansing and softening of encrusted dressing materials that are dried and adhered to the wound surface and difficult to remove. Actinol® Solution Plus softens the encrusted bandage and dressing, while eliminating the pathogenic microorganisms that have settled on them and prevents the wound from being infected.

Actinol® Solution Plus can be applied directly onto the wound as well as onto the bandage and gauzes. Materials

such as bandage, gauze, wiping clothes, etc. should not be used again. Actinol® Solution Plus can be heated to body temperature prior to use.

Usage with dressings and bandages Thoroughly wet dressings and bandages adhered to the wound surface with Actinol® Solution Plus, wait ca. 5-10 min. For softening and remove them carefully from wound. In order to maintain a moist wound climate, the dressing can be changed daily, if necessary, several times a day. Actinol® Solution Plus is compatible with all modern wound dressings.

Catheter care Possible infections can be avoided by regular cleaning catheter inlet ports, flexible connectors and its environment by Actinol® Solution Plus. Spray Actinol® Solution Plus on the area from a short distance and wait ca. 1-5 min. The ingredients of Actinol® Solution Plus are particularly mucous membrane-friendly, cleanse catheter and PEG material without damaging it. As with wound care, encrusted sticking plasters can be removed easily after wetting with Actinol® Solution Plus.

Side Effects Actinol® Solution Plus contains polyhexanide. Even though it is not irritating, it can cause allergic reactions in the bands and the patient's hands. Actinol® Solution Plus is reported that polyhexanide may cause an anaphylactic shock in rare cases (less than 1 out of 10,000). When allergic skin testing is recommended should not use. Patients hypersensitive to polyhexanide should not use the product. The biocompatibility index of Actinol® Solution Plus is found to be >1 and has a high tissue compatibility. Actinol® Solution Plus may cause allergic reactions such as itching (urticaria) and rash (eczema). In very rare cases, a mild burning sensation may occur after application of Actinol® Solution Plus, however this sensation disappears in a couple of minutes.

**Contraindications** Actinol® Solution Plus in the following cases: It is not suitable for use in the cartilage and aseptic joint surgeries, in the middle and internal ear and in the eye (Actinol® Solution Plus contacts the aseptic cartilage, wash immediately with Ringer solution or physiological saline solution).

- If the patient is known or suspected to be hypersensitive to any of the ingredients,
- On central nervous system or meningis,
- Along with anionic surfactants, cleansing soaps, creams, oils, enzymes, etc. These substances should be removed fully from the wound surface before use.
- For peritoneal lavage and washing,
- It should not be used in the first four months of pregnancy. Depending on the risk-benefit ratio, it should be used later in pregnancy should be evaluated according to the doctor's recommendation.

**Warnings and Precautions** Read the instructions carefully before use. Only for external application onto skin. Do not use for infusion or injection purposes. Do not swallow! Keep out of reach of children. Do not use if the sealing under the cap is torn or broken. Close the bottle immediately after use to avoid any contamination. The rim of the bottle should be protected from direct contact to avoid contamination while using. Bottles, which directly contact the wound or are contaminated in some other ways should be disposed. Keep only in original packaging contact from sunlight. Dispose according to local regulations.

**For newborns and infants:** Due to the insufficiency of clinical information, Actinol® Solution Plus should be used only under a strict and close medical supervision on newborns and infants.

This product is a Class II Medical Device according to the 93/42/EEC and the 2007/47/EC amendments.

**Ingredients:** Poloxamer, Polyhexamethylene biguanide (PHMB, 0.1%) and excipients.

**Appearance and odor:** Clear, colorless and almost odorless solution.

**Shelf life:** Please see the expiry date. Close the cap tightly and store in a cool, dry and frost-free environment. Shelf life is 2 (two) years at room temperature. Do not store an opened bottle above 35 °C. Consume the open product within 8 (eight) weeks.

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